

Sorvall ST 8 Compact Centrifuge Series

Fits in and stands out

Centrifuges

Fits in and stands out

Small benchtop centrifuges for clinical and research applications

The Thermo Scientific™ Sorvall™ ST 8 Compact Centrifuge Series is designed for clinical and research applications and provides exceptional capacity in a compact design. The centrifuges are available as ventilated or refrigerated models.



Fits in

- to your space, with its compact footprint. Ergonomic designs enhance placement flexibility in the lab.
- to your applications, accommodating multi-laboratory settings with 14 rotor options to adapt with evolving clinical and research needs, in both ventilated and refrigerated options.

Clinical applications

- Clinical chemistry
- Clinical microbiology
- Hematology
- Immunology

2



Research applications

- Cellular biology
- Microbiology
- Genomics / molecular biology
- Proteomics
- Biochemistry
- Pharmaceutical studies



Stands out

- designed for high safety, solid reliability and consistent results.
- featuring innovative Thermo Scientific™ technologies
- anticipating and addressing your most pressing lab needs even the ones you don't know you have yet.



Innovative technologies



Thermo Scientific™ Auto-Lock™ Rotor Exchange

Secure push-button application versatility and cleaning convenience



Thermo Scientific™ ClickSeal™ Biocontainment Lids

One-handed, certified¹ sample protection



Maximized swing-out capacity

Process up to 8×50 mL conical tubes, $24 \times 5/7$ mL blood tubes, microplates and microtubes, all in one versatile centrifuge



Fixed angle flexibility

Spin 50 mL conical tubes at over 12,000 x g or a wide range of microtubes at speeds up to 30,279 x g

Biocontainment certification by the Public Health Laboratory Service, Microbiology Services, Porton Down, UK or TÜV NORD CERT GmbH, Germany





Simple solutions Flexible processing

Secure, push-button **Auto-Lock** rotor exchange in as little as 3 seconds delivers:

- Trouble-free rotor installation and removal
 - No tools are required
 - No complicated rotor changing procedures
- · Clear chamber access for cleaning convenience
- Quick switch between applications from conical tubes to blood tubes, and from microtubes to microplates – and the flexibility to evolve with the changing needs of your laboratory

Optimize bench space with maximized capacity:

- 48 microtubes, 8 x 50 mL conical tubes or 4 microplates for research needs
- Micro-volume rotor portfolio supports 0.2, 0.25, 0.5, 1.5, 2.0,
 5.0 mL microtubes and spin columns
- Clinical processing up to 24 x 5/7 mL blood tubes per run

Reliable quality

Safe operation

- Biocontainment sealing options, including certified ClickSeal lids for glove-friendly, one-handed operation
- Manufactured with quality materials providing broad chemical resistance with exceptional strength for durability and reliability
- Low noise levels, including just 52 dBA at performance of 17,850 rpm/30,279 x g in refrigerated units, to ensure a peaceful laboratory
- Conforms to the latest clinical and safety standards, such as UL, CE and IVD²

Biocontainment certification by the Public Health Laboratory Service, Microbiology Services, Porton Down, UK or TÜV NORD CERT GmbH, Germany

² See full list of standards & certifications on page 14



Micro-volume rotor portfolio



Sample-ready

User-friendly

- One-touch operation with pre-saved protocols
- Highly visible backlit display for easy reading of parameters across the lab
- Optional indicators at end of run, including automatic lid opening, full flashing display and adjustable audible signal
- Glove- and detergent-friendly

 Multilingual instructions-English, Dutch, French, German, Italian, Russian and Spanish—on programming, run conditions, alerts and



Clinical applications

Outstanding capacity in a compact design with a smart and simple interface



Rotors for clinical applications

Thermo Scientific™ Rotor	Description	Stand-out features		mance
TX-150 Swinging Bucket Rotor	· ·		Ventilated	Refrigerated
	Designed for a broad range of high- speed, high-capacity swing-out processing needs, with a capacity of 24 x 5/7 mL blood tubes per run and flexibility with a wide range of adapters.	Certified¹ sample protection with ClickSeal biocontainment lids.	4,500 rpm 3,260 x g	4,500 rpm 3,260 x g
TX-100S/TX-100 Rotors				
	Dedicated for routine clinical samples processing.	Capacity of 16 tubes per run or spin 8 tubes with biocontainment sealing option.	4,500 rpm 3,260 x g	4,500 rpm 3,260 x g
HIGHConic III Rotor				
	High g-force fixed angle rotor for efficient sample separations.	Perform platelet poor plasma tests in as little as 1 minute.	8,700 rpm 10,155 x g	9,500 rpm 12,108 x g
8 x 50 mL Individually Sealed Fi	xed Angle Rotor			
	Designed for high-risk sample processing for a range of clinical and research tubes.	Individual sealed canisters for high-risk samples.	n/a	5,600 rpm 5,014 x g
CLINIConic Fixed Angle Rotor				
	Spin 30 blood tubes or 15 mL conicals in a fixed angle configuration with fast acceleration to optimize blood separation applications.	Auto-Lock rotor exchange for ultimate ease-of-use.	4,400 rpm 3,030 x g	4,400 rpm 3,030 x g
Hematocrit Rotor				
	Achieve red blood cell capillary separations.	Result reading chart provided.	13,300 rpm 16,810 x g	13,300 rpm 16,810 x g
M10 Swinging Bucket Rotor				
	Centrifuge 4 microplates or 2 midi- deepwell plates, ideally suited for quick spins. Spin up to 4 standard microplates with optional certified ¹ ClickSeal biocontainment lids.	Transportable microplate carriers, complemented by Auto-Lock rotor exchange.	4,400 rpm 2,576 x g	4,400 rpm 2,576 x g





Research applications

Innovative design supports multiple users and applications to stay one step ahead of your evolving research needs



- 8 x 50 mL cell culture tubes,
 4 microplates, 48 microtubes, or
 30 spin columns for research needs
- Auto-Lock rotor exchange for application flexibility — go from 50 mL tube to microplates to microtubes effortlessly





Rotors for research applications

Thermo Scientific	c rotor	Description	Stand-out features		rmance
			Jana Sat Idata Co	Ventilated	Refrigerated
TX-150 Swinging	Bucket Rote	Designed for a broad range of swing-out processing needs at high speed and high capacity with the flexibility of a wide range of adapters.	Dedicated bucket available for 8 conical tubes capacity or round buckets with certified¹ ClickSeal biocontainment lids.	4,500 rpm 3,260 x g	4,500 rpm 3,260 x g
M10 Swinging Bu	cket Rotor				
		Centrifuge 4 standard microplates or 2 deepwell plates. Spin up to 4 standard microplates with optional certified¹ ClickSeal biocontainment lids. Non-shielded design for quick access to microplates.	Transportable microplate carriers, complemented by Auto-Lock rotor exchange.	4,400 rpm 2,576 x g	4,400 rpm 2,576 x g
MT-12 Swinging E	Bucket Roto				
		Design capable of spinning up to 12 x 1.5/2 mL microtubes and filtration columns, ideal for small-volume research applications.	Individual swinging bucket rotors for centralized pellets within tubes.	13,000 rpm 16,438 x g	13,000 rpm 16,438 x g
HIGHConic III Fix	ed Angle Ro	otor			
*		Process 6 x 50 mL tubes in fixed angle configuration and at RCF above 10,000 x g in certified¹ sealed conditions.	45 degree fixed angle rotor ideal for cell pelleting using 50 mL or 15 mL conical tubes	8,700 rpm 10,155 x g	9,500 rpm 12,108 x g
MicroClick 18 x 5	Fixed Angle	e Rotor			
*		18-place capacity for 5.0 mL microtubes for high throughput and over 22,000 x g to reduce run times.	Lightweight and ergonomic ClickSeal biocontainment lid allows sample viewing.	n/a	14,000 rpm 22,351 x g
Micro-Volume Fix	ced Angle Ro	otors			
*	MicroClick 24 x 2	Capable of supporting micro- volume protocols, such as nucleic acid preparation, PCR reaction set-up and filtration	Run 24 x 2 mL samples at 30,279 x g.	16,000 rpm 24,328 x g	17,850 rpm 30,279 x g
*	MicroClick 30 x 2	columns, choose from the complete micro-volume rotor portfolio.	High-rim design of the MicroClick 30 x 2 rotor supports spin columns with open caps.	14,000 rpm 21,694 x g	14,000 rpm 21,694 x g
*	Microliter 48 x 2	_	Maximize throughput with 48-place microliter rotor.	11,800 rpm 15,256 x g	12,900 rpm 18,233 x g
*	8 x 8 PCR Strip	-	Process up to 8 PCR strips with certified¹ sample protection.	15,000 rpm 17,860 x g	15,000 rpm 17,860 x g

Extensive rotor portfolio

For application flexibility







Versatile swinging bucket rotors

TX-150 Rotor

M10 Microplate Rotor

MT-12 Rotor

	Rotor	Max tube	Vent	tilated	Refriç	gerated	
Description	capacity (places x volume, mL)	dimensions (Ø x L, mm)	Max Speed (rpm)	Max RCF (x g)	Max Speed (rpm)	Max RCF (x g)	Cat. No
TX-150 Swinging Bucket Rotor, 90°, R _{max} 144 mm (order buckets separately)	4 x 145	50 x 110	4,500	3,260	4,500	3,260	7500570
50 mL Conical buckets (unsealed, no adapter needed) (set of 4)	8 x 50	29.5 x 120					7500570
Round buckets (set of 4)	4 x 145	50 x 100					7500570
ClickSeal biocontainment lids for round buckets (set of 4)							7500570
Replacement O-Rings for lids (set of 4)							7500572
Adapters for 50 mL conical buckets (sets of 2)				,			
15 mL Conical tube	8 x 15	17 x 123					7500580
Adapters for Round Buckets (sets of 4)				'			
145 mL Bottle (75005734)	4 x 145	50 x 100					Direct fit
1 100 mL Round Bottom Open-Top Tube	4 x 100	45 x 117					7500573
1 50 mL round bottom tube ²	4 x 50	35 x 115					7500573
2 50 mL Conical or Skirted Tube	4 x 50	29.5 x 120					7500573
3 30 mL Sterilin™ Universal Container	4 x 30	25 x 120					7500574
4 15 mL Conical tube	8 x 15	17 x 122					7500573
4 11 mL IVF Tube	8 x 11	17 x 122					750057
5 13 mL Urine tube	16 x 13	17 x 110					750035
5 12 mL Blood collection tube (Greiner™)	16 x 12	17 x 110					750035
5 10 mL Blood Collection or 15 mL Corex™/Kimble™ Tube	16 x 15	17 x 110					750035
6 5/7mL mL Blood collection tube	24 x 5/7	13 x 110					750057
7 3/5 mL Blood collection tube	28 x 3/5	13 x 110					7500574
7 3/5 mL Cryotube	28 x 3/5	13 x 110					7500574
8 1.5/2 mL Microtube (or Microtainer™ tube)	40 x 2	11 x 65					7500574
Rotor Packages							
Cell Culture Package INCLUDES: TX-150 Rotor (75005701), Round Buckets (75005702), Adapters for 50 mL conical tubes (75005736)	4 x 50	29.5 x 120					7500576
High Capacity Cell Culture Package INCLUDES: TX-150 Rotor (75005701), Conical Buckets (75005703), Adapters for 15 mL conical tubes (75005808, 4 pks)	8 x 50	18 x 124					7500576
Clinical Rotor Package INCLUDES: TX-150 Rotor (75005701), Round Buckets (75005702), ClickSeal biocontainment lids (75005707), Adapters for blood collection tubes: 5/7 mL (75005739) and 10 mL (75003504)	24 x 5/7	18 x 124					7500576
M10 Microplate Rotor, 90°, R _{max} 119 mm order buckets separately)		n/a	4,400	2,576	4,400	2,576	7500570
Unsealed Buckets (set of 2)	4 Standard or	Height					7500572
Sealed Buckets (set of 2) R _{max} 109mm ❖	2 Deepwell	< 33 mm					7500572
MT-12 Swinging Bucket Rotor, 90°, R _{max} 87 mm	12 x 2	11 x 50	13,000	16,438	13,000	16,438	750056
Replacement ClickSeal Biocontainment Lid (each)							7500573
Replacement Sealings (set of 2 with grease)							7500572

 $^{^{2}\,\}mbox{Used}$ together with adapter 75007542

























TX-100 Rotor

Clinical swinging bucket rotors

	Rotor capacity	Max tube	Venti	ated	Refrigerated		
Description	(places x volume, mL)	dimensions (Ø x L, mm)	Max Speed (rpm)	Max RCF (x g)	Max Speed (rpm)	Max RCF (x g)	Cat. No.
TX-100S Clinical Swinging Bucket Rotor with Sealed Carriers, 90°, R _{max} 144 mm	8 x 15	17 x 106	4,500	3,260	4,500	3,260	75005704
TX-100 Clinical Swinging Bucket Rotor with Carriers, 90°, R _{max} 144 mm	16 x 15	17 x 121	4,500	3,260	4,500	3,260	75005705
Adapters for TX-100/TX-100S Clinical Rotors (e	ach)						
13 mL Urine tube (TX-100 only)	16/8 x 13	17 x 110					Direct fit
10mL Blood collection tube	16/8 x 10	16 x 100					Direct fit
Sterilin Urine tube (TX-100 only)	16/8 x 13	16 x 110					75008817
Sarstedt Blood tube (TX-100 only)	16/8 x 10	16 x 75					75008818
5/7 mL BD Hemogard™/BD Vacutainer™ Tube	16/8 x 5/7	13 x 110					11172596
5 mL BD Hemogard tube	16/8 x 5	13 x 75					11172595
3 mL Blood collection tube	16/8 x 3	11 x 70					11172287
1.5/2 mL Microtube (or Microtainer™ tube)	16/8 x 1.5/2	10 x 41					11172288











HIGHConic III Rotor

8 x 50 mL Rotor

CLINIConic Rotor

MicroClick 18 x 5 Rotor

Fixed angle rotors

	Rotor	Max Tube	Vent	tilated	Refrig	gerated	
Description	Capacity (places x volume, mL)	Dimensions (Ø x L, mm)	Max Speed (rpm)	Max RCF (x g)	Max Speed (rpm)	Max RCF (x g)	Cat. No.
HIGHConic III Rotor, 45°, R _{max} 120 mm *	6 x 50	30 x 121	8,700	10,155	9,500	12,108	75005709
Replacement lid (each)							75005731
Replacement sealings (set of 2 with grease)							75003058
Adapters for HIGHConic III Rotor (set of 2)							
50 mL Round Bottom Tube and	6 x 50	30 x 115					Direct fit
50 mL Conical Bottom Tube 38 mL Round bottom tube	6 x 38	25.5 x 110					75005802
16 mL Round bottom tube	6 x 16	18 x 123					75005802
15 mL Conical tube	6 x 15	17 x 123					75005803
12 mL Round bottom tube	6 x 12	16 x 95					75005808
6.5 mL Round bottom tube	6 x 6.5	13.5 x 114					75005805
5 mL Conical microtube	6 x 5	17 x 100					75005805
3.5 mL Round bottom tube	12 x 3.5	11 x 100					75005776
1.5/2 mL Microtube	12 x 2	11 x 40					75005807
8 x 50 mL Individually Sealed Rotor 25°, R _{max} 143 mm	8 x 50 Sealed	30 x 121	n/a	n/a	5,600	5,014	75003694
Replacement biocontainment lids (set of 2)							75003011
Replacement O-ring kit							75003789
Adapters for 8 x 50 mL Individually Sealed	Rotor (set of 2)		I			,	
50 mL Round bottom tube	6 x 50	30 x 115					Direct fit
38 mL Round bottom tube	6 x 38	25.5 x 110					75005802
16 mL Round bottom tube	6 x 16	18 x 123					75005803
15 mL Conical tube	6 x 15	17 x 123					75005808
12 mL Round bottom tube	6 x 12	16 x 95					75005804
6.5 mL Round bottom tube	6 x 6.5	13.5 x 114					75005805
5 mL Conical microtube	6 x 5	17 x 100					75005770
3.5 mL Round bottom tube	12 x 3.5	11 x 100					75005806
1.5/2 mL Microtube	12 x 2	11 x 40					75005807
CLINIConic Rotor, 37°, R _{max} 140 mm	30 x 15	16.5 x 131	4,400	3,030	4,400	3,030	75003623
Adapters for CLINIConic Rotor (each)							
15 mL Round/conical bottom tube	30 x 15	16.5 x 131					Direct fit
10 mL Round bottom tube	30 x 10	16.5 x 95					75008817
5/7 mL BD Hemogard/ BD Vacutainer tube	30 x 5/7	13 x 106					11172596
5 mL BD Hemogard tube	30 x 5	13 x 75					11172595
MicroClick 18 x 5 Rotor, R _{max} 102 mm	18 x 5	17 x 62	n/a	n/a	14,000	22,351	75005765
Sealing set							75005726
Rubber seal grease							76003500
Adapters for MicroClick 18 x 5 (set of 2)							
1.5/2 mL Microtube	18 x 1.5/2	11 x 45					75005756







MicroClick 30 x 2 Rotor



MicroClick 24 x 2 Rotor



8 x 8 PCR Strip Rotor



Hematocrit Rotor

Fixed angle rotors

	Rotor Capacity	Max Tube	Vent	tilated	Refrig	jerated	
Description	(places x volume, mL)	Dimensions (Ø x L, mm)	Max Speed (rpm)	Max RCF (x g)	Max Speed (rpm)	Max RCF (x g)	Cat. No.
Microliter 48 x 2 Rotor, 45°, R _{max} 98 mm *	48 x 1.5/2	11 x 45	11,800	15,256	12,900	18,233	75003602
Replacement biocontainment lid (each)							70904727
Replacement O-Ring for lid (each)							75003349
MicroClick 30 x 2 Rotor, 45°, R _{max} 99 mm *	30 x 1.5/2	11 x 45	14,000	21,694	14,000	21,694	75005719
Replacement ClickSeal biocontainment lid (each)							75005730
Replacement O-Ring for lid (each)							75005726
MicroClick 24 x 2 Rotor, 45°, R _{max} 85 mm *	24 x 1.5/2	11 x 45	16,000	24,328	17,850	30,279	75005715
Replacement ClickSeal biocontainment lid (each)							75005725
Replacement O-Ring for lid (each)							75003405
Adapters for Microliter 48 x 2, MicroClick 30 x	2 and MicroClick	24 x 2 Rotors (s	ets of 30)				
0.2 mL PCR tube	48/30/24 x 0.2	6.5 x 20					75005752
0.5 mL Microtube	48/30/24 x 0.5	8 x 44					75005753
0.25 mL Microtube	48/30/24 x 0.25	6 x 46					75005754
8 x 8 PCR Strip Rotor, 45°, R _{max} 71 mm	8 x 8 PCR strips	n/a	15,000	17,860	15,000	17,860	75005720
Replacement ClickSeal biocontainment lid (each)							75005730
Replacement O-Ring for lid (each)							75005726
Hematocrit Rotor, 90°, R _{max} 85 mm, includes linear scale	24 Capillaries	n/a	13,300	16,810	13,300	16,810	75005733
Capillaries (pack of 100)							76000923
Sealing kit							75000964
Replacement linear scale reading graph							76000938



Sorvall ST 8 Centrifuges







Product specifications

Specificatio	ns	Sorvall ST 8	Sorvall ST 8R	Sorvall ST 8FR				
Swing Bucket	Max Capacity	4 x 145 mL (with TX-150 rotor)	4 x 145 mL (w	vith TX-150 rotor)				
Rotors	Max Speed	4,500 rpm (with TX-150 rotor)	4,500 rpm (w	ith TX-150 rotor)				
	Min Speed	300 rpm (with TX-150 rotor)	300 rpm (wit	th TX-150 rotor)				
Max RCF		3,260 x g (with TX-150 rotor)	3,260 x g (with TX-150 rotor)					
Fixed Angle Max Capacity		6 x 50 mL (with HIGHConic rotor)	6 x 50 mL (with HIGHConic rotor)					
Rotors	Max Speed	16,000 rpm (with MicroClick 24 rotor)	17,850 rpm (with	MicroClick 24 rotor)				
	Min Speed	300 rpm (with MicroClick 24 rotor)	300 rpm (with MicroClick 24 rotor)					
	Max RCF	24,328 x g (with MicroClick 24 rotor)	30,279 x g (with	MicroClick 24 rotor)				
Control System	n		Microprocessor					
Drive System			Direct, brushless induction drive					
Rotor Locking	System		Auto-Lock rotor exchange					
Imbalance Det	ection System	Continuous	vibration measurement, with rotor ma	ss correction				
Programs		4 direct program buttons	3 direct program buttons, 96 additional programs accessible via folder					
Pulse (Short) F	lun	Yes						
Acceleration/D	eceleration Rates	2/2 (standard and soft)						
Centrifuge Cha	amber	Brushed stainless steel						
Timer Range			10 second to 99h 59min, continuous					
Temperature R	lange	n/a	-10°C to 40°C					
Pre-Cooling Fu	ınction	n/a	Yes, with direct button					
Refrigeration S	System	n/a	CFC free					
Sound Level		<58 dBA (with TX-150 rotor) <61 dBA (with MicroClick 24 rotor)	<56 dBA with all rotors <52 dBA (with MicroClick 24 rotor)					
Other Features	S		Certified1 ClickSeal biocontainment lids					
		Multilingual selection	: English, Dutch, French, German, Italian	, Russian and Spanish				
Max Power Co	nsumption	310 W	750 W	850 W				
Dimensions	Height (lid open)	67 cm (26.8 in)	70 cm (27.5 in)	117-122 cm (45.6-47.6 in), adjustable				
	Height (lid closed)	31 cm (12.2 in)	32 cm (12.6 in)	79-84 cm (30.7-32.7 in), adjustable				
	Width	37 cm (14.6 in)	46 cm (18.1 in)	47 cm (18.3 in)				
	Depth	48 cm (18.9 in)	67 cm (26.4 in)	52 cm (20.5 in)				
Weight		35 kg (77 lb)	71 kg (156 lb)	87 kg (191 lb)				
Standards &	Standards	IEC 61010-	IEC 61010-1, IEC 61010-2-020, IEC 61010-2-101, EN 61326-1					
Certifications	Certifications	UL Listed, C	E marked, IVD compliant, Certified bio	ocontainment ¹				

Ordering information

Product	Cat. No.
Sorvall ST 8 centrifuge, ventilated, 120 V, 60 Hz	75007200
Sorvall ST 8 centrifuge, ventilated, 230 V, 50/60 Hz	75007201
Sorvall ST 8 centrifuge, ventilated, 100 V, 50/60 Hz	75007202
Sorvall ST 8R centrifuge, refrigerated, 120 V, 60 Hz	75007203
Sorvall ST 8R centrifuge, refrigerated, 230 V, 50/60 Hz	75007204
Sorvall ST 8R centrifuge, refrigerated, 100 V, 50/60 Hz	75007205
Sorvall ST 8FR floor-standing centrifuge, refrigerated, 120 V, 60 Hz	75007207
Sorvall ST 8FR floor-standing centrifuge, refrigerated, 230 V, 50/60 Hz	75007208
Sorvall ST 8FR floor-standing centrifuge, refrigerated, 100 V, 50/60 Hz	75007206

Complete packages

For your application needs

Application package	Accessories included	Centrifuge	Cat. No.
Multi-discipline			
	TX-150 rotor TX-150 buckets (set of 4) TX-150 ClickSeal biocontainment lids (set of 4)	Sorvall ST 8R refrigerated centrifuge, 120 V	75230906
	 Adapters 4 x 50 mL conical Adapters 8 x 15 mL conical HIGHConic III 6 x 50 rotor 	Sorvall ST 8 ventilated centrifuge, 120 V	75200906
	M10 Microplate rotorMicroClick 24 x 2 rotor	Sorvall ST 8R refrigerated centrifuge, 230 V	50144456
High capacity cell culture			,
Second .	• TX-150 rotor - TX-150 oval buckets (set of 4)	Sorvall ST 8R refrigerated centrifuge, 120 V	75230395
	- Adapters 8 x 15 mL conical	Sorvall ST 8 ventilated centrifuge, 120 V	75200395
		Sorvall ST 8R refrigerated centrifuge, 230 V	50147009
		Sorvall ST 8 ventilated centrifuge, 230 V	50147005
Cell culture			
	TX-150 rotor TX-150 buckets (set of 4) TX-150 ClickSeal biocontainment lids (set of 4) Adapters 4 x 50 mL conical Adapters 8 x 15 mL conical	Sorvall ST 8R refrigerated centrifuge, 120 V	75233637
		Sorvall ST 8 ventilated centrifuge, 120 V	75203637
		Sorvall ST 8R refrigerated centrifuge, 230 V	50144426
	- Adapters of X To THE conteat	Sorvall ST 8 ventilated centrifuge, 230 V	50140771
Routine clinical			
	 TX-150 rotor TX-150 buckets (set of 4) TX-150 ClickSeal biocontainment lids (set of 4) Adapters 24 x 5/7 mL blood 	Sorvall ST 8R refrigerated centrifuge, 120 V	75230102
		Sorvall ST 8 ventilated centrifuge, 120 V	75200102
	collection - Adapters 16 x 10 mL blood	Sorvall ST 8R refrigerated centrifuge, 230 V	50144430
	collection	Sorvall ST 8 ventilated centrifuge, 230 V	50140786
PPP separation			'
	HIGHConic III 6 x 50 rotor Adapters for 5/7 mL	Sorvall ST 8 ventilated centrifuge, 120 V	75200947
Molecular biology			·
	• TX-150 rotor - TX-150 buckets (set of 4)	Sorvall ST 8R refrigerated centrifuge, 120 V	75230115
	- TX-150 ClickSeal biocontainment lids (set of 4)	Sorvall ST 8 ventilated centrifuge, 120 V	75200115
	 Adapters 4 x 50 mL conical Adapters 8 x 15 mL conical MicroClick 24 x 2 rotor 	Sorvall ST 8R refrigerated centrifuge, 230 V	50147011
		Sorvall ST 8 ventilated centrifuge, 230 V	50147008







Learn more at thermofisher.com/benchtopcentrifuges

For Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific uses or applications. © 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Greiner is a registered trademark of Greiner. Corex is a registered trademark of Corning Inc. Kimble is a registered trademark of Kimble Glass. BD Hemogard and BD Vacutainer are registered trademarks of BD Biosciences. EXT6231 1223

thermo scientific