# VersaCool

## Refrigerated Bath Circulator

Manual Part Number U01265, April 4, 2016

**Multilingual Setup Wizards**

**Multilingual Essential Safety Instructions**

**Installation**

**Operation**

**Preventive Maintenance**

**Troubleshooting**

# Table of Contents

|  |            |
|--|------------|
| <b>Preface</b> .....                       | <b>P-1</b> |
| Compliance .....                           | P-1        |
| Unpacking .....                            | P-2        |
| After-Sale Support .....                   | P-2        |
| Sample Nameplate.....                      | P-2        |
| Feedback.....                              | P-3        |
| <b>Chapter 1 Safety</b> .....              | <b>1-1</b> |
| Safety Factors.....                        | 1-1        |
| Personal Protective Equipment.....         | 1-3        |
| Training .....                             | 1-3        |
| <b>Setup Wizards</b> .....                 | <b>W-1</b> |
| <b>Chapter 2 General Information</b> ..... | <b>2-1</b> |
| Description .....                          | 2-1        |
| Intended Use .....                         | 2-1        |
| Specifications .....                       | 2-2        |
| <b>Chapter 3 Installation</b> .....        | <b>3-1</b> |
| Bath Installation .....                    | 3-1        |
| Ventilation.....                           | 3-1        |
| Electrical Requirements.....               | 3-2        |
| Connectivity .....                         | 3-3        |
| Lid Reversal/Removal .....                 | 3-4        |
| External Circulation .....                 | 3-5        |
| Approved Fluids .....                      | 3-6        |
| Additional Fluid Information .....         | 3-8        |
| Tubing Requirements .....                  | 3-10       |
| Optional Auto Refill.....                  | 3-11       |
| Nitrogen Purge.....                        | 3-11       |
| Filling Requirements.....                  | 3-12       |
| Draining .....                             | 3-12       |
| <b>Chapter 4 Operation</b> .....           | <b>4-1</b> |
| VersaCool Touchscreen .....                | 4-1        |
| Initial Start Up .....                     | 4-2        |
| Daily Start Up.....                        | 4-2        |
| Changing a Set Point .....                 | 4-4        |
| Shut Down.....                             | 4-4        |
| Touchscreen Displays.....                  | 4-5        |
| Viewing/Changing a Settings Screen.....    | 4-6        |
| Viewing/Changing a Home Screen .....       | 4-10       |
| Viewing/Changing a Service Screen.....     | 4-18       |
| High Temperature Cutout.....               | 4-21       |

|                   |  |            |
|-------------------|--|------------|
| <b>Chapter 5</b>  | <b>Preventive Maintenance.....</b>             | <b>5-1</b> |
|                   | Condenser Fins.....                            | 5-1        |
|                   | Hoses.....                                     | 5-1        |
|                   | Electrical Power Cord .....                    | 5-1        |
|                   | Cleaning.....                                  | 5-2        |
| <b>Chapter 6</b>  | <b>Troubleshooting .....</b>                   | <b>6-1</b> |
|                   | Error Messages.....                            | 6-1        |
|                   | Checklist .....                                | 6-3        |
|                   | Testing the Safety Features .....              | 6-4        |
| <b>Appendix A</b> | <b>AC Serial Communications Protocol.....</b>  | <b>A-1</b> |
| <b>Appendix B</b> | <b>Analog IO.....</b>                          | <b>B-1</b> |
| <b>Appendix C</b> | <b>Mobile Communications.....</b>              | <b>C-1</b> |
| <b>Appendix D</b> | <b>GNU LESSER GENERAL PUBLIC LICENSE .....</b> | <b>D-1</b> |

**Declaration of Conformity**

**Warranty**

# Preface

## Compliance

### US/Canada

UL Listed to US and Canadian Standards (File E473332)

### European Union - CE

The Declaration of Conformity is available on request.

### WEEE

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with 'wheelie bin' symbol:

Thermo Fisher Scientific has contracted with one or more recycling/ disposal companies in each EU Member State, dispose of or recycle this product through them. Further information on Thermo Fisher Scientific's compliance with these Directives is available at:

[www.thermoscientific.com/WEEERoHS](http://www.thermoscientific.com/WEEERoHS)

## Unpacking

The bath is supplied with an electrical power cord. Do not discard the packaging until the cord is located and the bath is operating.

If the bath shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

**The bath does not have handles. Take into account its weight, 80 pounds (36 kilograms), when unpacking and transporting. We recommend two people lift it from the bottom.**

**Position baths in an upright position for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor.**

## After-sale Support


Thermo Fisher Scientific is committed to customer service both during and after the sale. If you have questions concerning the bath operation, or questions concerning spare parts or Service Contracts, call our Sales, Service and Customer Support phone number, see next page for contact information.

Before calling, please obtain the following information:

- bath BOM number
- bath serial number
- bath firmware version (see Chapter 4)
- power source voltage

The bath's BOM and serial number are located on the nameplate label on the rear of the bath, see label sample on the next page.

## Sample Nameplate

## Feedback

We appreciate any feedback you can give us on this manual. Please e-mail us at:

[tcmanuals@thermofisher.com](mailto:tcmanuals@thermofisher.com)

Please include the manual part number and the revision date listed on the front cover.

### **Thermo Fisher Scientific**

25 Nimble Hill Road  
Newington, NH 03801  
Tel : (800) 258-0830 or  
(603) 436-9444  
Fax : (603) 436-8411  
[www.thermoscientific.com/tc](http://www.thermoscientific.com/tc)

### **Sales, Service, and Customer Support**

25 Nimble Hill Road  
Newington, NH 03801  
Tel: (800) 258-0830  
Sales: 8:00 am to 5:00 pm  
Service and Support: 8:00 am to 6:00 pm  
Monday through Friday (Eastern Time)  
Fax: (603) 436-8411  
[service.tc.us@thermofisher.com](mailto:service.tc.us@thermofisher.com)

# Chapter 1 Safety

## Safety Factors

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your bath. If you have any questions concerning the operation of your bath or the information in this manual, please contact us. See inside cover for contact information.

**DANGER** indicates an imminently hazardous situation that, if not avoided, will result in death or serious injury.

**WARNING** indicates a potentially hazardous situation that, if not avoided, could result in death or serious injury.

**CAUTION** indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the bath's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates the presence of hot surfaces.

This label indicates read the manual.

**Observe all warning labels.**

**Never remove warning labels.**

**The bath's construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided.**

**The circuit protector located on the rear of the bath is not intended to act as a disconnecting means.**

**The bath's power cord is used as the disconnecting device, it must be easily accessible at all times.**

**Never operate the bath with damaged cords.**

**The bath is not intended for use with Personal Protection Ground Fault Interrupter (GFI) outlets with a rating of 10 mAmp or below.**

**If GFI outlets are required at the point of installation, Equipment Protection GFI outlets with a rating above 10 mAmp are recommended.**

**Never place the bath in a location or atmosphere where excessive heat, moisture, or corrosive materials are present.**

**Many refrigerants are undetectable by human senses and are heavier than air replacing the oxygen in an enclosed area causing loss of consciousness. Refer to the R134A refrigerant's SDS for additional information.**

**Never connect the process fluid inlet or outlet fittings to the building water supply or any water pressure source.**

**Never operate the bath without fluid in the reservoir.**

**Use only the approved fluids listed in this manual.**

**To prevent freezing, never operate the bath below 5°C with only water in the reservoir.**

**Other than water, before using any approved fluid, or when performing maintenance where contact with the fluid is likely, refer to the manufacturer's SDS and EC Safety Data sheet for handling precautions.**

**Ensure the fluid will not generate toxic gases. Flammable gases can build up over the fluid during usage.**

**Do not use automotive antifreeze. Commercial antifreeze contains silicates, or any Organic Acid Technologies (OATs), that can damage the pump seals. Use of automotive antifreeze voids the manufacturer's warranty.**

**Never use corrosive fluids with this bath. Use of these fluids also voids the manufacturer's warranty.**

**When using a process fluid mixture of ethylene glycol and water or propylene glycol and water, check the fluid concentration and pH on a regular basis. Changes in concentration and pH can impact system performance. See Chapter 3.**

**Never operate damaged or leaking equipment.**

**Never operate the bath or add fluid to the reservoir with panels removed.**

**Do not clean the bath with solvents, only use a soft cloth and water.**

**Always drain the bath before moving. Drain the bath before it is transported and/or stored. Store the bath in the temperature range -25°C to 60°C (with packaging), and <80% relative humidity.**

**Always turn the bath off and disconnect the supply voltage from its power source before moving the bath or before performing any service or maintenance procedures. Transport the bath with care. Sudden jolts or drops can damage its components.**

**Refer service and repairs to a qualified technician.**

**Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and voids the manufacturer's warranty and safety compliance.**

## **Personal Protective Equipment**

There are no special personal protective equipment requirements needed to perform normal operation unless required by the manufacturer of the fluid you are using. We do recommend wearing eye protection and gloves.

## **Training**

The user must review and understand all the sections in this manual before operating the bath.

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Pekskärmen visar kort:

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# Setup Wizard

This Setup Wizard is intended for initial start up only. For all other procedures, or if these steps are not clear, you must refer to the manual.

Read Safety Factors in Chapter 1 before starting.

### Approved Fluids:

- Filtered/single distilled water (pH 7-8)
- Deionized water (Maximum 1 MΩ-cm, compensated)
- Distilled water with Nalco biocide and inhibitor
- Distilled water with chlorine (5 ppm)
- 50/50 Laboratory Grade EG/Water
- 50/50 Laboratory Grade PG/Water
- SIL 180      SIL 200

### What you need to get started:

- Adjustable wrenches
- Appropriate hose or plumbing
- Appropriate size hose clamps
- Approved fluid

Leave baths in an upright position at room temperature for 24 hours before starting. This ensures the lubrication oil drains back into the compressor.

- The plumbing connections for external circulation are located on the rear of the bath.  is the return flow from the external application.  is the outlet flow to the external application (supply side). The connections are male M16 x 1. Remove the union nuts and plates to install the supplied 1/4", 1/2", 8 mm or 12 mm hose barbs. Also supplied are 1/4" MNPT and 1/2" MNPT hose barbs used primarily with quick disconnects.

To prevent damage to the bath's plumbing, use a 19 mm backing wrench when removing/installing the external connections.

- Ensure the reservoir drain port is closed and that the process discharge and return fittings are plumbed or capped. Ensure any residue is thoroughly removed from the reservoir before filling.

- To avoid spilling, place your samples/trays into the bath before filling.

- Slowly fill the reservoir to between the **MIN** and the water or oil **MAX** fill lines.

- After the power cord is connected place the circuit protector located on the rear of the bath to the **I** position.


- The touchscreen momentarily displays:



- Then the Setup Wizard Display appears.


- Touch Setup Wizard to display Language. Touch the up/down arrow until the desired setting appears and then touch  to save the selection and to sequence through the other touchscreens or  to return to the previous screen.

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

Name your bath There is a 20 character limit. The name identifies the bath when downloading files.

- After several seconds the Power touchscreen appears.

- If desired, touch  to start the bath or refer to the manual to adjust the bath's settings.
-  is replaced by . Touch  when you are ready to stop the bath.
- After starting check all plumbing connections for leaks.
- Adjust the High Temperature Cutout (HTC) safety device, refer to the manual.

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## Chapter 2 General Information

### Description

All Thermo Scientific VersaCool™ Refrigerated Bath Circulators can pump to an external system. All have a digital display and easy-to-use touchscreen, five programmable set point temperatures, acoustic and optical alarms, and offer adjustable high temperature protection.

### Intended Use

Intended use is to provide heating, cooling or constant temperature to applications placed into the reservoir and/or pumped to an external application in accordance with all the procedures and requirements stated in this manual.

Intended for use in a laboratory environment, on a bench top only.

Typical lab applications include:

|                            |              |
|----------------------------|--------------|
| Analytical Instrumentation | Lasers       |
| General Laboratory Cooling | Condensers   |
| Rotary Evaporators         | Bio-reactors |
| Microscopes                | Histology    |


The serial tag on the rear identifies the bath, its electrical requirements and the refrigeration data.

## Specifications

| VersaCool   |  |
|---|--|
| <b>Process Fluid Temperature</b> °C<br><b>Range</b> °F                        | -20 to 150<br>-4 to 302  |
| <b>Ambient Temperature</b> °C<br><b>Range</b> °F                              | 10 to 40<br>50 to 104  |
| <b>Internal Bath Stability</b> °C<br>Water<br>EG/Water<br>PG/Water            | ±0.03<br>±0.06<br>±0.05  |
| <b>Cooling Capacity</b> watts   | 425  |
| <b>Heating Capacity</b> watts<br>115V<br>230V                                 | 1200<br>2000   |
| <b>Bath Volume</b> liters<br>gallons  | 7.0<br>1.8   |
| <b>Bath Work Area Dimensions</b><br>(D x W x L) cm<br>inches                  | 15.5 x 17.3 x 29.9<br>6.0 x 6.8 x 11.7   |
| <b>Net Weight</b> kg<br>lb  | 36<br>80   |
| <b>Pumping Capacity</b><br>Max flow rate (lpm/gpm)<br>Max pressure (mbar/psi) | 14.5/3.8<br>1200/17.40   |
| <b>Electrical Requirements</b> (VAC/Hz)<br>(Voltage ±10%)                     | 100-120 /50 or 60<br>or<br>200-240 /50 or 60   |
| <b>Interfaces</b>   | Remote Sensor Port<br>Micro USB<br>USB<br>RS232 RS485<br>Bluetooth<br>Ethernet<br>Smart-View<br>Multifunction I/O Port |
| <b>Compliance</b>   | CE<br>RoHS<br>cULus<br>WEEE  |
| <b>Maximum Relative Humidity</b> °C<br>(Non Condensing)°F                     | 80% for temperatures up to 31°(88°F)<br>decreasing linearly to 50% relative humidity at 40°C (104°F)                   |
| <b>Operating Altitude</b> meters<br>feet                                      | Sea Level to 2000<br>Sea Level to 6560   |
| <b>Overvoltage Category</b>   | II   |
| <b>Pollution Degree</b>   | 2  |
| <b>Degree of Protection</b>   | IP20   |
| <b>R134A Refrigerant</b> kilograms<br>ounces                                  | 0.18<br>6.35   |
| <b>Storage Temperature Range</b> °C<br>°F                                     | -25 to +60<br>-13 to +140  |
| <b>Location of Use</b>  | Indoor only  |

Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C to +24°C (unless otherwise specified) ambient condition, at nominal operating voltage. Other fluids, process or ambient temperatures, altitude or operating voltage will affect performance. Pump specifications are nominal values of ±10%. Reduce work area depth to 10.3 cm (4 inches) when using silicon oil.


- Thermo Fisher Scientific reserves the right to change specifications without notice.



Specifications at sea level using at a +20°C ambient condition, at nominal operating voltage. Other fluids, process or ambient temperatures, altitude or operating voltage will affect performance. The bath was operating with the lid closed and no heat load lid or any external circulation.

Minimum temperatures are only achieved with the above conditions. If your conditions vary from above, minimum temperature will most likely not be achieved.


- Thermo Fisher Scientific reserves the right to change specifications without notice.



Time to Temperature, Water


Time to Temperature, Water


Time to Temperature, EG/Water 50/50


Time to Temperature, EG/Water 50/50


Time to Temperature, PG/Water 50/50


Time to Temperature, PG/Water 50/50



Time to Temperature, Sil 180


Time to Temperature, Sil 200


centimeters  
inches

Shipping crate dimensions:  
 H = 78.11cm / 30.75"  
 W = 38.74cm / 15.25"  
 D = 78.11cm / 30.75"  
 Shipping weight = 46.3kg / 102lbs

# 3

## Chapter 3 Installation

### Bath Installation

The bath is designed for continuous operation and for indoor use.

**Never place the bath in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present.**

**Baths should be left in an upright position for 24 hours at room temperature before starting. This will ensure the lubrication oil has drained back into the compressor.**


### Ventilation

Inadequate ventilation will result in excessive internal bath temperatures that may result in system shutdown and component failures.

Operating with water up to 90°C and up to 40°C ambient, or Silicon Oil up to 150° and up to 35°C ambient, requires a minimum clearance of 0 inches on one side, open on the other side and 6 inches on the rear.

**Operating with Silicon Oil above 125°C and 35°C ambient or above requires a minimum clearance on one side of 12", open on the other side and 6" in the rear.**

Minimum Clearance

A = 6" (15.2 cm)

B = 0"

B = 12" (30.5 cm) for Sil Oil above 125°C and 35°C ambient or above

C = Must be open

## Electrical Requirements

Refer to the nameplate on the rear of the bath for specific electrical requirements.

The bath is intended for use on a dedicated outlet.

**The bath construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided.**

**EN 61000-3-11 Compliance:** This equipment complies with EN/IEC 61000-3-11 provided that main power system impedance ( $Z_{max}$ ) is less than or equal to  $0.10\Omega + j0.06\Omega$  at the interface point between the user's supply connection and the public system. Consult with the local supply authority if necessary for system  $Z_{max}$  determination.

The circuit protector on the rear of the bath is designed to protect the bath's internal components.

**Note** If the circuit protector activates allow the temperature to cool before resetting. Restart the bath. Contact us if the circuit protector activates again.

**The bath's electrical power cord is used as the disconnecting device, it must be easily accessible at all times.**

**Ensure the cord does not come in contact with any of the plumbing connections, reservoir contents or tubing.**

**Note** Before inserting the electrical cord into the connection ensure the circuit protector is in the **0** (off) position.

Once the cord is connected to the bath, connect the other end to the main power source.


## Connectivity

**Never apply line voltage to any of these connections.**

This connection located on the rear of the bath is used to accommodate a temperature monitoring probe including the optional Smart-View system probe.

### Remote Temperature Sensor

The remote temperature sensor on the rear of the bath requires a 4-pin connector that must mate to a LEMO # ECP.1S.304.CLL. The bath uses a 3 wire sensor, but a 4 wire sensor can be used (pins 3 and 4 are interconnected in the control head). The pin-out is:

Pin 1 and 2 = Pt100 +                      Pin 3 and 4 = Pt100 -

See Chapter 4 for instructions to enable the remote sensor.


### Ethernet

Indicates the location of the standard RJ45 Ethernet connection for use with future communication/control protocols.

### Serial Communications

The RS232 9-pin connector on the front and rear of the bath are used to select and activate serial communications. There is also a RS485 9-pin connector on the rear of the bath for use with future communication/control protocols. See Appendix A for additional information.


#### RS232

| Pin #     | Function            |
|-----------|---------------------|
| 1, 4, 6-9 | No connection       |
| 2         | TX                  |
| 3         | RX                  |
| 5         | GND = Signal ground |


TX = Transmitted data from bath

RX = Received data to bath

#### RS485

| Pin # | Function      |
|-------|---------------|
| 1-7   | No connection |
| 8     | T+            |
| 9     | T-            |

### Multifunction Port

The 15-pin Multifunction port on the rear of the bath has contacts for fault/warnings detection, remote on/off and can set and report temperature to your PLC. See Appendix B for additional information.

### Auto Refill

Indicates the location of the connection for the optional auto refill accessory.

In addition there are USB, RS232 and MicroUSB connections located on the front of the bath.

The USB connection is used for data upload, e.g., firmware upgrades, user ramping programs, and data download, e.g., data logging and user defined ramping programs. It will be used with future communication/control protocols.

The RS232 and MicroUSB connections will also be used with future communication/control protocols.


### Mobile Communications


The VersaCool Mobile App can use any iOS or Android smart phone or tablet to control and monitor up to seven VersaCool baths. Mobile communications is enabled when the circuit protector on the rear of the bath is in the **I** position. See Appendix C for additional information.

### Lid Reversal/Removal

To remove or reverse the lid slide the tab on the side of the lid to the unlocked position and then lift the lid. Once reversed and properly positioned use the tab to lock the lid into place.

**Locked**

**Unlocked**

## External Circulation

The plumbing connections for external circulation are located on the rear of the bath.

 is the return flow from the external application.  is the outlet flow to the external application (supply side). The connections are male M16 x 1. Remove the union nuts and plates to install the supplied 1/4", 1/2", 8 mm or 12 mm hose barbs. Also supplied are 1/4" MNPT and 1/2" MNPT hose barbs used primarily with quick disconnects.

**Note** When not using external circulation the plumbing connections must be capped.

**Note** To prevent damage to the bath's plumbing, use a 19 mm backing wrench when removing/installing the external connections.

## Approved Fluids

**Only use the approved fluids listed in this manual.**

Thermo Fisher Scientific takes no responsibility for damages caused by the selection of an unapproved fluid.

**Ensure, when selecting the heat transfer fluid, that no toxic gases are generated. Inflammable gases can build up over the fluid during usage. Other than water, handle and dispose of approved fluids in accordance with the fluid manufacturer's specification and/or the fluid's Safety Data Sheet (SDS).**

**To set the fluid temperature alarms, use the bath's menu selection to identify the fluid, see Chapter 4.**

**When using water above 80°C closely monitor the fluid level, frequent top-offs will be required. It will also create steam.**

**Glycol/Water mixtures require top-offs with pure water, otherwise the percentage of glycol will increase resulting in high viscosity and poor performance.**

**Ensure the over temperature cut-off point is set lower than the fire point for the heat transfer fluid selected.**

**The highest working temperature, as defined by the EN 61010 (IEC 1010), must be limited to 25°C below the fire point of the fluid.**

Refer to Chapter 4 for fluid set point ranges.

### **Water, Distilled Water or Deionized Water** (1 MΩ-cm maximum)

Normal tap water leads to calcareous deposits necessitating frequent bath decalcification.

Calcium tends to deposit itself on the heating element. The heating capacity is reduced and service life shortened.

Well or municipal water should not be used unless it has been tested by a lab to meet the minimum requirements shown on the next page. The primary risks of using untested and untreated tap water are corrosion caused by an out of range pH value and/or the concentration of dissolved solids, and those dissolved solids precipitating out of solution and forming scale.

| Water Quality and Standards                         |   |                 |
|---|---|-----------------|
| Process Fluid                                       | Permissible (PPM)                           | Desirable (PPM) |
| <b>Microbiologicals</b><br>(algae, bacteria, fungi) |   |                 |
|   | 0   | 0               |
| <b>Inorganic Chemicals</b>                          |   |                 |
| Calcium   | <25   | <0.6            |
| Chloride  | <25   | <100            |
| Copper  | <13   | <10             |
|   | 0.020 ppm if fluid in contact with aluminum |                 |
| Iron  | <0.3  | <0.1            |
| Lead  | <0.015                                      | 0               |
| Magnesium   | <12   | <0.1            |
| Manganese   | <0.05                                       | <0.03           |
| Nitrates\Nitrites                                   | <10 as N                                    | 0               |
| Potassium   | <20   | <0.3            |
| Silicate  | <25   | <1.0            |
| Sodium  | <20   | <0.3            |
| Sulfate   | <25   | <1              |
| Hardness  | <17   | <0.05           |
| Total Dissolved Solids                              | <50   | <10             |
| <b>Other Parameters</b>                             |   |                 |
| pH  | 6.5-8.5                                     | 7-8             |
| Resistivity   | 0.01*                                       | 0.05-0.1*       |

\* MΩ-cm (compensated to 25°C)

Unfavorably high total ionized solids (TIS) can accelerate the rate of galvanic corrosion. These contaminants can function as electrolytes which increase the potential for galvanic cell corrosion and lead to localized corrosion such as pitting. Eventually, the pitting may become so extensive that refrigerant will leak into the reservoir.

As an example, raw water in the United States averages 171 ppm (of NaCl). The recommended level for use in a water system is between 0.5 to 5.0 ppm (of NaCl).

Recommendation: Initially fill the reservoir with distilled or deionized water. Do not use untreated tap water as the total ionized solids level may be too high. This reduces the electrolytic potential of the water and prevent or reduce the galvanic corrosion observed.

| Advantages                     | Disadvantages  |
|--------------------------------|--|
| Readily available              | Can support bio-growth   |
| Great capacity to carry heat   | Evaporates readily at top of temperature range   |
| Low viscosity                  | Corrosive to ferrous metals  |
| Compatible with many materials | pH needs to be maintained  |
| Distilled water is sterile     | Distilled and DI water need to be tested   |
| -                              | Distilled and DI water must be purchased or require specialized equipment to produce on site |
| -                              | Deionizing does not guarantee sterility  |

## Chlorine

Short term usage of tap water may not cause any adverse affects on the bath or your application, but in the long term problems may arise. To help alleviate these problems Thermo Fisher Scientific recommends the use of chlorine.

The duration of time that chlorine remains in solution depends on factors such as water temperature, pH and availability of direct sunlight. We recommend maintaining chlorine levels at proper levels using chlorine test strips, generally 1 to 5 ppm is adequate.

For best results, maintain the pH of the fluid between 6.5 and 7.5. Do not add additional chlorine without first determining the concentration ratio that already exists in the fluid supply. Corrosion and degradation of the circulation components can result from concentration ratios that are too high. Contact our customer support for additional information.

## Nalco

We recommend Nalco algaecide/corrosion inhibitor for applications with aluminum and/or where biological growth is a concern.

## 50/50 Laboratory Grade Ethylene Glycol with Water, by volume 50/50 Laboratory Grade Propylene Glycol with Water, by volume

Glycols are the most common freeze-point suppressant. **Note** Low-end temperatures using 50/50 EG/Water or 50/50 PG/Water were achieved using lower pump speeds.

Excess glycol deteriorates the temperature accuracy due to its high viscosity.

| Advantages  | Disadvantages   |
|---|---|
| Wider temperature range than water alone                              | Does not remove heat as well as water alone   |
| Low cost and readily available  | The viscosity goes up as the temperature goes down (PG is more viscous than EG at lower temperatures) |
| Good capacity to carry heat   | Can support bio-growth  |
| Low viscosity   | Evaporates readily at top of temperature range  |
| Compatible with many materials  | Can become corrosive without inhibitor  |
| Sterile (PG is Less toxic/safer to use, easier to dispose of than EG) | pH needs to be maintained to prevent corrosion  |
| -   | Toxic - needs to be handled and disposed of with care   |

## Silicon Oil (Sil 180, Sil 200)

Silicone oil comes in a variety of viscosities and temperature ranges. The higher-viscosity silicone oils typically also have a higher fire point and work well at temperatures that cannot be achieved with water (> 95°C). The lower viscosity silicone oils are intended for use below 0°C, where their lower viscosity allows them to still be effectively pumped or circulated. Silicone oil can also be used between 5 °C and 95 °C for applications that cannot use water.

| Advantages  | Disadvantages                                    |
|---|--|
| Wide temperature range                            | Expensive  |
| Low rate of evaporation                           | Messy – fumes can create oil film                |
| No bio-growth                                     | Disposal   |
| Non-corrosive                                     | Can become corrosive without inhibitor           |
| Widely compatible                                 | Not compatible with silicone hose or ABS plastic |
| Long life span (over 1 year) and negligible smell | Immiscible in water (does not mix)               |

## Additional Fluid Information

- Other than water, do not use any approved fluid until you have read and understood the manufacturer's instruction for use and the Safety Data Sheet (SDS).
- Ensure any fluid residue or any other material is thoroughly removed before filling the reservoir with a different fluid.
- Always wear protective clothing, especially a face shield and gloves.
- Avoid spattering on any of the circulating bath's components, always slowly add fluid. When adding, point the opening of a container away from yourself.
- For proper ventilation, use a fume hood.
- Do not allow any ignition sources in the vicinity.
- When using oil:

Oil contaminated with water may effervesce when oil is heated above 100°C causing the reservoir to overflow. In order to dry the oil, operate the bath at 100°C for one hour and then at 105°C for another hour.

Clean and dry the reservoir and covers, especially gable covers. Any item placed in the reservoir must also be dry.

Vacuum out reservoir, drain line and process lines.

## Tubing Requirements

**Ensure none of the tubing comes in contact with the power cord.**

Tubing is normally used to connect the pump to an external application.

**Note** The maximum allowable length of tube depends largely on the size, form and material of the external reservoir. The length of tube and its diameter, combined with the circulating capacity, have a large affect on the temperature stability. Whenever possible, use a wider tube diameter and place the application as close as possible to the bath.

**Operating the bath at either the high or low end temperature limits will lead to similar temperatures on the tube surface, this is even more critical with metal tubing.**

- The required tube material depends on the heat transfer liquid used
- Tubes must not be folded or bent
- After prolonged use, tubes may become brittle or they may get very soft, check them on a daily basis and replace if necessary
- Secure all tube connections using clamps

When using the internal reservoir only, the plumbing connections must be closed with the supplied plate and union nuts.

**Ensure the tubing you select meets your maximum temperature and pressure requirements.**

### Plastic and rubber tubing

If plastic and rubber tubes are used, ensure that the tubes selected are fully suitable for the particular application, i.e., that they will not split, crack or become disengaged from their connections.

Connect the tubing using the supplied tube fittings. They are attached to the plumbing connections with a supplied coupling nut.

Supplied Fittings:

M16 - to - 1/4" NPT (M)

8mm Hose Barb

1/4" Hose Barb

M16 - to - 1/2" NPT (M)

12mm Hose Barb

1/2" Hose Barb

We highly recommend using foam rubber insulation on the tubing and the fittings.

## Metal tubing

Thermo Scientific metal tubing (stainless steel insulated) offers a particularly high degree of safety and is suitable for both low and high temperatures/liquids.

The metal tubing is attached directly to the plumbing connections, gaskets are not required.

**Note** Do not subject tubing to mechanical strain and ensure any specified bend radius is not exceeded.

Tubing is available in lengths of 0.5, 1.0 and 1.5 meters. Couplings for connecting tubes are also available.

The smallest opening inside the metal tubes is 10 mm. The metal tubing is provided with coupling nuts (M16 x 1, DIN 12 879, part 2) at either end.

## Optional Auto Refill Accessory

The auto refill accessory attaches to the  connections on the rear of the bath. The auto refill provides make-up fluid to replace bath fluid lost to evaporation, etc. It will require a pressurized fluid source.

When make-up fluid is available and the bath is turned on, for normal operation the auto refill will keep the reservoir fluid at the levels set using the bath's touchscreen, see Chapter 4. The auto refill shuts off when the fluid reaches those levels, or if the bath is turned off.

**Note** The maximum input pressure rating is 150 psi.

## Nitrogen Purge

Baths are equipped with nitrogen purge line and vent designed to accept a constant flow of dry nitrogen into the reservoir. The nitrogen blankets the cooling fluid reducing fluid evaporation.

Connect your nitrogen line to the **N2** 1/8" OD tube on the rear of the bath.

**Note** Limit the pressure to no more than 3 PSI.

## Filling Requirements

Ensure the reservoir drain port on the front of the bath is closed and that all plumbing connections are securely plumbed or capped. Also ensure any residue is thoroughly removed from the reservoir before filling.

**Before using any fluid refer to the manufacturer's MSDS and EC safety data sheets for handling precautions.**

To avoid spilling, consider fluid displacement when objects are placed in the reservoir.

*Slowly* fill the reservoir until the fluid is between the MIN and the water or oil MAX fill lines.

**Note** Operating the bath with the fluid below the MAX fill line will reduce cooling capacity.

**Avoid overfilling, oil-based fluids expand when heated. The maximum fill when using Silicon Oil is 6 liters (1.65 gallons).**

When pumping to an external system, keep extra fluid on hand to maintain the proper level in the circulating lines and the external system.

**Note** Monitor the fluid level whenever heating the fluid.

## Draining

**Always drain the bath before moving or storing.**

**Before draining any fluid refer to the manufacturer's MSDS and EC safety data sheets for handling precautions.**

**Ensure the fluid is at a safe handling temperature, ~40 °C. Wear protective clothing.**

- Remove the lower front panel to expose the drain.
- If desired, attach an 8 mm id tube on the drain.
- Place a suitable receptacle underneath the drain.
- Slowly turn the drain plug until flow is observed.

**Note** When the reservoir empties there is still 1/3 liter of fluid in the bath lines. To ensure all the fluid is removed ensure the flow has totally stopped from the drain plug before closing.

- Use a wet-vac on the drain connection to remove 100% of the process fluid.
- When the flow stops close the drain plug.

# 4

## Chapter 4 Operation

### Initial Start Up

For initial start up refer to the Setup Wizard in the front portion of this manual.

### Daily Start Up

**Before starting, double check all connectivity, electrical and plumbing connections.**

Do not run the bath until fluid is added to the reservoir. Have extra fluid on hand. If the bath does not start refer to Chapter 6 Troubleshooting.

- Place the circuit protector located on the rear of the bath to the I position.
- The touchscreen momentarily displays:

**Note** The bath may “click” as it loads.

When loading is complete the **Power** touchscreen will appear.

 This icon indicates the bath is using the internal temperature sensor. The procedure to switch to the remote,  , or outlet sensor,  , is addressed later in this chapter.

 This icon indicates the approximate reservoir fluid level.

• Touch  to place the bath into the sleep mode, this causes the screen to go blank. Touch any part of the blank screen to come out of the sleep mode.

The display also shows the process fluid temperature as well as the set point. The set point is the desired process fluid temperature.

• Touch  to start the bath. The  icon becomes  , touch it to stop the bath.

Additional icons can be displayed but they are on only when that component is running/enabled.

 This icon indicates the pump is running, it is always on when the bath is running.

 This icon indicates the heater is on, it flashes when the heater is cycling on and off.

 This icon indicates the refrigeration is on.

 This icon indicates the timer, used to start/stop the bath at a specific time, is enabled.

 This icon indicates a connectivity option is enabled.

Touch one of the four icons on the bottom of the screen to operate and configure the bath.

**Note** The bottom boarder of the selected icon is highlighted with a green glow.

**Power** is used to bring up the screen needed to start and shut down the bath. It also displays the process fluid temperature, set point (program name if a program is running), temperature sensor and reservoir fluid level.

**Settings** is primarily used to configure the bath's displays, e.g., date, time and language.

**Home** is primarily used to configure the bath's operating values, e.g., set points, alarms, fluid type and pump speed.

**Service** is primarily used to maintain/troubleshoot the bath, e.g., perform a calibration, display error messages and reset factory values.

The screen can also display a keypad and keyboard.

The keypad automatically appears when a numerical value needs changing. Touch the keys to enter the desired setting and then touch **OK** to accept the value and close the keypad. Touching **X** before touching **OK** also closes the keypad but ignores the change.

Invalid entries turn red when **OK** is touched.

The keyboard automatically appears to enter/edit text. Touch the keys to enter the desired text and then touch **Enter** to accept the change and close the keyboard. Touching **X** before touching **Enter** also closes the keyboard but ignores the change.

Touch here for capital letters  
Touch here for numerals

When keypad and keyboard are used  and  always appear at the bottom of the display. You must touch either  to save the change, or touch  to not save the change, in order to return to the previous screen.

Except for Full Screen Charting, if any screen is not touched for 60 seconds a status screen appears. This screen displays the selected process fluid and the process fluid temperature. Touching the status screen returns you to the previous display.

## Changing the Set Point

Touch **28.00°C** to bring up the set point selection display.

Either touch **XX.X°C** to bring up the keypad to make changes or

touch  next to the desired preset set point to select it.

Touch  or  to return to the **Power** touchscreen.

For set point ranges and how to change the five preset set points refer to Viewing/Changing a Home Screen in this chapter.

# Touchscreen Displays

Touching any one of the four icons at the bottom of the screen brings up access to a lower level of displays used to operate the bath. Some of the displays have additional lower levels.



**Power**

**Settings**

**Home**


**Service**


## Viewing/Changing a Settings Screen

Viewing/changing a Settings screen requires you touch  or  to exit the screen. Save the change by touching . To not save the change touch .

Use the **Date** screen to view/change the date.

If a screen has a down arrow, , touch it to remove the four icons.

If desired, touch the up arrow, , to display the four icons.

Use the **Time** screen to view/change the time and format.

Use the **Display** screen to adjust the screen intensity. Touch and drag  to the desired intensity. Touch  to bring up the keypad to enter the exact desired value.

The factory default is **Max**.

**Note** Low-end temperatures using 50/50 EG/Water or 50/50 PG/Water were achieved using lower pump speeds.

Use the **Units** screen to change the displayed temperature scale, °C or °F. The selected scale is highlighted in green.


The factory default is **°C**.

Use the **Language** screen to change the language displayed on the touchscreen.


- Chinese
- English
- French
- German
- Italian
- Japanese
- Spanish

The factory default is **English**.

Use the **Resolution** screen to change the displayed temperature resolution, 0.1 or 0.01. The selected resolution is highlighted in green.

The factory default is **0.1**.

Use the **User Access** screen to open or lock user access. **Open** allows full access to all functionality.

**Locked** allows partial access to only select a program or a pre-determined set point. You cannot change alarms, warnings or faults.

The factory default is **Open**.

To create a password and lock the bath touch **Set/Reset** and use the keypad to create a five digit password and then touch **OK**. Reenter the password to confirm and then touch **OK** again.

To open the bath touch **Open** and enter the password. To re-lock the bath touch **Locked** and enter the password. If you forget your password call us with the VersaCool's serial number.

Use the **Timer** screen to set and enable/disable the timer. Use it to set the desired time for the bath to turn on. The screen also displays the current time.

When enabled the On indication is green, when disabled the Off indication is red.

The factory default is **Off**.

Use the **Auto Restart** screen to enable/disable auto restart.

When auto restart is enabled and the bath shuts down as a result of a power failure, when power is restored the bath automatically restarts and operates at the saved values. **Consider any possible risks before enabling.**

The factory default is **Off**.

Use the **Sleep Mode** screen to enable/disable the sleep mode. Once enabled, if the Status Screen is showing and no icon is touched for 60 seconds the screen will enter Sleep Mode and powers down. Touching the screen in Sleep Mode returns it to the Status Screen.

The default is **Off**.

Use the **PID** screen to view/adjust the bath's heat PID values.

Thermo Fisher recommends that only a qualified technician adjust the PID values. Incorrect values will hamper bath performance.

Factory preset values for water:  
**P** = 5.0, **I** = 0.30 and **D** = 0.00


Use the **Name** screen to name your bath. There is a 20 character limit.

The name identifies the bath when downloading files.

The default is **VersaCool 1**.

## Viewing/Changing a Home Screen

Making a change requires you to save the change by touching . To return to the previous screen and not save the change touch .

Use the **Set Points** screen to verify/adjust the set points and RTA. The set point is the desired process fluid temperature. Touch the desired SP icon and then touch its value window to bring up the keypad.

The factory default values are:

**SP1** = 20°C, **SP2** = 40°C,  
**SP3** = 50°C, **SP4** = 70°C  
**SP5** = 95°C.

You cannot adjust the set point closer than 0.1° to either of the fluid's system limits shown below.

The fluid is selected using the **Fluid Type** screen, see next page.

| Set Point Temperature Range* |                        |                                   |
|------------------------------|------------------------|-----------------------------------|
| Approved Fluid               | Altitude (meters/feet) | Range                             |
| Water                        | <610/2000              | +5°C to +90°C (+41°F to +194°F)   |
| Water                        | >610/2000              | +5°C to +85°C (+41°F to +185°F)   |
| Water                        | >1370/4500             | +5°C to +80°C (+41°F to +176°F)   |
| PG/Water                     | <1070/3500             | -16°C to +90°C (3°F to +194°F)    |
| PG/Water                     | >1070/3500             | -16°C to +85°C (3°F to +185°F)    |
| EG/Water                     | All                    | -18°C to +80°C (0°F to +176°F)    |
| SiI 200                      | All                    | +20°C to +135°C (+68°F to +275°F) |
| SiI 180                      | All                    | -10°C to +150°C (14°F to +302°F)  |


\*Maximum process fluid temperatures are limited based on altitude.

If the displayed temperature does not accurately reflect the actual temperature in the reservoir a Real Time Adjustment (RTA) is required. The RTA can be set  $\pm 10^{\circ}\text{C}$  ( $\pm 18^{\circ}\text{F}$ ).

Each set point has its own RTA.

As an example, if the temperature is stabilized and displaying 20°C but a calibrated reference thermometer reads 20.5°C, set the RTA to -0.5°C. After you enter an RTA value allow the display to stabilize before verifying the bath temperature. If display accuracy is required, we recommend repeating this procedure at various set point temperatures and on a regular basis.

Use the **Alarms** screen to set the temperature warning and fault limits and to enable/disable the temperature and fluid level audible alarms.

The speaker icon, , is visible when the alarm is enabled. The mute speaker icon, , is visible when disabled.


Alarm default settings are determined by the fluid type, see **Fluid Type** screen.

High Temp Fault = +3 above maximum fluid temperature

High Temp Warning = +3 maximum fluid temperature

Low Temp Warning = -3 below minimum fluid temperature

Low Temp Fault = -3 below minimum fluid temperature

Use the **Fluid Type** screen to select the desired fluid. The screen also displays each fluid's temperature limits.


The factory default is **Water**.

If you change fluid type and the current set point is outside one of the new fluid type temperature limits, a warning will appear asking **“Do you want to change Fluid Type?”**

Touching **Yes** brings up the current set point screen with a keypad to input the new set point.

Touching **No** returns the **Fluid Type** screen.

**Note** With Sil 200, heater power is reduced to 60% when connected to 230V.

Use the **Pump Speed** screen to set the desired pump speed percentage, 40% - 100%. Slide the button along the bar or touch the button to type in an exact speed percentage using the keypad.


**Note** When Sil 180 or Sil 200 is the selected fluid the pump speed range is limited to 90% - 100%.

Use the **Charting** screen to view a graphical display of the bath's performance.

Touch the temperature sensor and/or **Set Point** window to display them. Touch **+** and **-** to adjust the displayed temperature range and/or time interval.

**Note** The current time is displayed on the right side of the display.

**Note** When running a program the set point value changes in order to meet the program's step time to temperature parameters.

Touch the **Full Screen** icon  to enlarge the screen.

Touch  to return to the above display.

Use the **Program** screen to **CREATE** a program or **DELETE, VIEW** or **RUN** a previously created program. (There may be a screen delay while the programs load.)

Touch the desired program's name to delete, view or run it.

Define your program as a series of set points with a known period of time between each. Each interval is one step of the program. Pay careful attention to the first part of your program. What conditions must exist at the beginning of your process? For example, at the starting set point you may wish to program an initial period of constant temperature to allow for thermal stabilization.

Also consider the bath's limitations when designing programs. Temperature or time parameters which exceed the bath's performance capabilities will result in unsatisfactory operation. If reaching the ramp set point temperatures is important, you will have to operate the bath between the desired set points and note the duration before programming the ramp.

It is possible to create a program calling for very rapid changes in temperature. Although the bath may not be capable of producing such changes, it may be practical to program such steps as a way to cause the fastest possible temperature change.

Touch **CREATE** to start building a new program. There can be up to 10 Programs, each with up to 30 Steps.

Touch **NEW** to define a **STEP**.

For **STEP 1** touch the vacant window below **START TEMP** and enter the desired value and then touch the vacant window below **END TEMP**. By default, the end temp for any step becomes the start temp for the following step. Changing the start temp for the following step will result in an error message when the program is saved.

Next enter the **TIME**. The maximum duration for each **STEP** is 999 minutes (~16.5 hours).

**# OF LOOPS** and **LOOP TO STEP** enables program step repetition. For example, after STEP 4 set the # OF LOOPS to 3 and LOOP TO STEP to 1. Once the program initially reaches the end of STEP 4, it will loop back to STEP 1 (LOOP TO STEP) and run through STEPS 1 through 4 again. In this instance the program will loop back three times. The entire sequence of steps will run a total of four times. After the fourth time the program will sequence to STEP 5. The maximum number of loops is 100. Enter 0 in both fields if looping is not required.

**Note** When using this feature ensure the loop end temperature agrees with the start temperature. In this example, ensure the end temperature for Step 4 is the same as the start temperature for Step 1 or you will receive a step error message when trying to save the program.

Touch **NEW** to create additional steps. Touch a **STEP** row to edit/delete it.

**Variance** is used to set a starting temperature range, the program starts when the fluid temperature is within this range. For example, if the desired **Start Temp** is 25°C and the **Variance** is set to 0.5°C, the program automatically starts only when the temperature is between 24.5°C to 25.5°C.

The program has an optional **Soak** feature. When enabled, this feature pauses the program timer until the temperature reaches set point,  $\pm$  variance. This assures the temperature reaches set point for all the steps before the program continues to the next step.


Touch **SAVE** to save the program into the bath's non-volatile memory.

Touch **RUN** to start the program and bring up the Program Display.

**Note** The program *will not start* until the process fluid temperature is at the **START TEMP**  $\pm$  the **Variance**.

Touch **Full Screen** to enlarge the display.

Use the **Connectivity** screen to select the desired option. The respective icon will highlight.

Mobile communications is always activated, see Appendix C.

**Note** You cannot save changes using the touchscreen when Analog IO is enabled.

Touch the desired option in order to configure its settings.

**Baud Rate:** 300, 600, 1200, 2400, 4800, 9600, 19200, 38400

**Parity:** Odd or Even

**Data Bits:** 7 or 8

**Stop Bits:** 1 or 2


**Addresses:** 1 to 99 (RS 485 only)

Refer to Appendix A for serial communication protocol information.

Refer to Appendix B for details on Analog IO.

Use the **Data Log** screen to view/adjust the data logging settings.

The factory default is **Off**.

**On Off** = User wants to start/stop Data Logging

**Sample Rate** = Time between data points (in seconds)

**Collection Period** = Total time the user wants Data Logging to run (in minutes)


Example:

Sample Rate = 30 seconds

Collection Period = 60 minutes

When the Data Logging is complete, the user will have 121 Data Points.


(121 includes the Data Point at “Time = Zero”)


Use the **Data Transfer** screen to upload/download data files between the bath and a storage device.

Uploading/downloading includes data recording, error logs, charting data and ramp data.

Touchscreen firmware upgrades can also be uploaded.


Use the **Control Probe** screen to select which temperature probe monitors/controls the bath. **Internal** is located in the bath, **Outlet** is located on the external circulation fluid out port. **Remote** is the sensor at the application.

The factory default is **Internal**.

Selecting **Remote** allows you to set fluid temperature limits **Above** and **Below** the set point.

With limit on the VersaCool will only operate within the set range regardless of what is happening at the external application. This is a safety function and protects the bath from trying to achieve undesirable temperatures.

The factory default is **5°C**.

Use the **Auto Refill** screen to enable and configure the optional auto refill.

Touch **Time Out (minutes)** to set the maximum time the option will operate. Setting the time to 0 disables the option.

**Note** Turn the auto refill off before adding or removing objects into the bath's reservoir.

Touch **Level Selection** for the auto refill range of operation.


**Note** The auto refill uses the selected fluid to determine the **Full** level. **Empty** does not let the reservoir run dry, there will still be fluid in the reservoir.

## Viewing/Changing a Service Screen

Making a change requires you to save the change by touching . To return to the previous screen and not save the change touch .

Use the **Configuration** screen to view the bath's configuration. This information cannot be changed.

Use the **Calibration** screen to calibrate the **Internal**, **Outlet** and/or **Remote** temperature sensor.

This procedure requires a calibrated reference thermometer.

Before starting ensure the selected sensor is controlling the temperature and that its RTAs is set to 0.

Touch the desired sensor.

**Note** To speed up the **Internal** calibration procedure drain the reservoir to the **MIN** fill level. For accuracy place your reference thermometer at the rear of the reservoir, as close to the screen as possible.

To calibrate the **Internal** or **Remote** sensor:

- Touch the **Sp x** value window to bring up the keypad and enter either the desired high or desired low calibration temperature. Keep in mind the selected fluid temperature limits.
- Allow the **Current** temperature to stabilize near the calibration temperature.
- Once stabilized, for three to five minutes, touch either **High** or **Low** value window and enter the temperature displayed on your reference thermometer and then touch **Cal** on the right of **High** or **Low**.
- Touch **Sp x** again and enter the other calibration temperature, repeat the procedure.
- Once both points are entered, touch  to complete the calibration.

Contact us if you need to calibrate the **Outlet** temperature sensor.

Use the **Run Times** screen to view component total running time. This information cannot be changed.

Use the **Diagnostics** screen to view component status. This information cannot be changed.

Use the **Error Messages** screen to view component error messages. This information cannot be changed although you can clear all the messages by touching **Clear Screen**.

Use the **Factory Reset** screen to return the five stored set points, all the alarms, the PID values and the pump speed back to their factory settings.

Touch **RESET** and then **Confirm** or **Cancel** the procedure.

This screen appears once the reset is complete.

## High Temperature Cutout

To protect your application, the adjustable High Temperature Cutout (HTC) ensures the heater will not exceed temperatures that can cause serious damage. A High Temperature fault occurs when the temperature of the sensor exceeds the set temperature limit.

In the event of a fault the bath shuts down and displays a fault message, see Chapter 6. Identify and correct the cause of the fault must before restarting the bath. A primary reason for the HTC to trip is a low reservoir fluid level.

The HTC is factory preset fully clockwise to the highest possible setting. To set the cutout start the circulator and adjust the set point a few degrees higher than the highest desired fluid temperature. Allow the circulator to stabilize at the temperature set point. Then, using a flathead screwdriver, slowly turn the red dial counterclockwise until the circulator shuts down and the fault message appears.

Before you can restart the bath it has to cool down a few degrees.

To restart the bath press the black reset ring surrounding the red dial. If Auto Restart is enabled the bath restarts, if disabled use the Daily Start Up procedure.

**Note** We recommend periodically rechecking operation or if the bath is moved.

# 5

## Chapter 5 Preventive Maintenance

**Laboratory Grade Ethylene Glycol (EG) is poisonous and flammable. Before performing any preventive maintenance refer to the manufacturer's most current SDS and EC Safety Data sheet for handling precautions.**

**Disconnect the power cord prior to performing any maintenance.**

**Handle the bath with care, sudden jolts or drops can damage its components.**

There are no user serviceable components within the bath's panels.

### Condenser Fins

In order to maintain cooling capacity, clean the condenser fins two to four times per year, depending on the operating environment.

Switch off the bath circulator and unplug the power cord.

Clean the fins with compressed air.

### Hoses

Inspect and tighten the external hoses and clamps daily.

### Electrical Power Cord

Ensure any replacement cord is properly rated.

## Cleaning

**Wear protective clothing and take appropriate measures when handling the cleaning agent.**

Before cleaning the bath's surfaces, to protect labels, the nameplate, electrical connections, painted and plastic surfaces and to prevent the cleaning agent from entering through any vent openings, mask off all areas except the reservoir.

Clean the bath's surface only with a soft cloth and warm water.

After time, the circulating bath's stainless steel surfaces may show spots and become tarnished. Normal stainless steel cleaners can be used.

Clean the reservoir and built-in components at least every time the bath liquid is changed. Use only water and a soft cloth.

**Do not use scouring powder or any substance containing solvents.**

The inside of the bath circulator must be kept clean in order to ensure a long service life. Quickly remove substances containing acidic or alkaline substances and metal shavings as they could harm the surfaces causing corrosion. If corrosion (e.g., small rust marks) occur in spite of this, cleaning with stainless steel caustic agents has proved to be suitable. Apply these substances according to the manufacturer's recommendations.

To clean the touchscreen use a microfiber cloth. For dirtier screens add a small amount of plain water (soap is not necessary).

# 6

## Chapter 6 Troubleshooting

### Error Messages

An error message indicates an unusual condition. WARNING and FAULT error messages are a result of exceeding one of the bath's Alarm limit settings (see Chapter 4), exceeding a sensor factory preset safety value or a safety switch is activated. With either message the alarm, if enabled, will sound.

In the case of a WARNING the bath, if running, will continue to run. Touch  to see if the message clears and, if enabled, silence the alarm. If the limit was only temporarily exceeded the message will not reappear.

In the case of a FAULT the bath will shut down but the screen will continue to flash the message. Touch  to clear the screen and, if enabled, silence the alarm. Once the cause of the shut down is identified and corrected, restart the bath. If the cause of the fault was not corrected the message will reappear.

| Error Messages            |  |
|---------------------------|--|
| Message                   | Cause/Action   |
| <b>HTC Fault</b>          | <p>high temp protection limit exceeded</p> <p>turn red dial full clockwise</p> <p>allow bath to cool down</p> <p>reset HTC by pressing black ring</p> <p>restart the bath</p> <p>reset HTC to desired setting, see Chapter 4</p> <p>when operating at high temperatures ensure pump is in a high-speed range</p> <p>if the HTC fault cannot be cleared, the bath must be serviced by an authorized Thermo Scientific Temperature Control Service Technician.</p>  |
| <b>High Temp Fault</b>    | <p>adjustable high temp fault protection limit exceeded</p> <p>check limit setting</p> <p>check fluid selection</p> <p>ensure bath has adequate ventilation</p>  |
| <b>High Temp Warning</b>  | <p>adjustable high temp warning protection limit exceeded</p> <p>check limit setting</p> <p>check fluid selection</p> <p>ensure bath has adequate ventilation</p>  |
| <b>Low Temp Warning</b>   | <p>adjustable low temp warning protection limit exceeded</p> <p>check limit setting</p> <p>check fluid selection</p>   |
| <b>Low Temp Fault</b>     | <p>adjustable low temp fault protection limit exceeded</p> <p>check limit setting</p> <p>check fluid selection</p>   |
| <b>High Fluid Level</b>   | <p>high level protection limit exceeded</p> <p>check fluid level</p> <p>drain fluid, if required</p>   |
| <b>Low Fluid Level</b>    | <p>low level protection limit exceeded</p> <p>check fluid level</p> <p>check for leaks</p>   |
| <b>Internal RTD Open</b>  | open internal temp sensor  |
| <b>Internal RTD Short</b> | internal temp sensor short   |
| <b>Remote RTD Open</b>    | open remote temp sensor  |
| <b>Remote RTD Short</b>   | remote temp sensor short   |
| <b>RA RTD Open</b>        | open refrigeration temp sensor   |
| <b>RA RTD Short</b>       | refrigeration temp sensor short  |
| <b>RA RTD Temp Fault</b>  | <p>check voltage supply</p> <p>the refrigeration may need servicing</p>  |
| <b>HPC Fault</b>          | <p>high pressure protection limit exceeded</p> <p>check for obstructions to air flow</p> <p>the refrigeration may need servicing</p>   |

| Known HMI Software Issues (Version 3D)  |  |
|---|--|
| Issue   | Action   |
| <b>Unable to update internal and external RTA values</b>                        | Manually adjust set point until fluid temperature matches the reference thermometer.           |
| <b>Temperature alarm values do not update when fluid type is changed</b>        | Wait 15 seconds after changing fluid before accessing the Alarm screen.                        |
| <b>Pump speed goes to 100% whenever fluid type is changed.</b>                  | Manually reset pump speed.   |
| <b>Programs start immediately (doesn't wait to get to the start temp value)</b> | Reprogram the start of the program to include time to get to the desired starting temperature. |
| <b>Long programs do no complete</b>   | Split into smaller programs.   |
| <b>HMI locks up if error message displayed for greater than three minutes</b>   | Recycle circuit protector on the rear of the bath.   |
| <b>Analog-IO calibration does not work</b>                                      | None   |

## Checklist

### Bath will not start or shuts down

- Check display for error codes, see Error Displays in this chapter.
- Ensure stop wasn't accidentally pressed.
- Ensure the circuit protector is in the on ( I ) position.
- Check the line cord connection to your power supply and at the bath.
- Make sure supply voltage is connected and matches the bath's nameplate rating  $\pm 10\%$ .
- Restart the bath.

### No display

- Touch the screen, it may be in the Sleep Mode.
- Cycle the bath's circuit protector.

### Temperature display reads 320°C

- Remote temperature probe selected but there is no probe attached to the bath.
- Attach a remote probe or select the internal temperature sensor.

### All set points display 0°C

- Cycle the bath's circuit protector.

continued on next page

### Bath will not circulate process fluid

- Check the reservoir level. Fill, if necessary.
- Check the application for restrictions in the cooling lines.
- The pump motor overloaded. The pump's internal overtemperature overcurrent device shuts off the pump causing the flow to stop. Possible causes are low fluid, debris in system, operating bath in a high ambient temperature condition or excessively confined space.
- Allow time for the motor to cool down.
- Make sure supply voltage matches the bath's nameplate rating  $\pm 10\%$ .

### Inadequate temperature control

- Verify the set point.
- Make sure the condenser is free of dust and debris.
- Check the fluid concentration.
- Ensure bath installation complies with the site requirements in Chapter 3.
- Make sure supply voltage matches bath's nameplate rating  $\pm 10\%$ .
- If the temperature continues to rise, make sure your application's heat load does not exceed the rated specifications.
- Check for high thermal gradients (e.g., the application load turning on and off or rapidly changing).
- Low-end temperatures using 50/50 EG/Water or 50/50 PG/Water were achieved using lower pump speeds.

### Program will not advance

- If **Soak** is enabled the program *will not advance* until the process fluid temperature is at the **END TEMP**  $\pm$  the **Variance**.

### Unable to save a changes using the touchscreen


- Check the start up screen to see if the communication icon is on,  .
- If it is on ensure Analog IO is not enabled,  . You cannot save a changes using the touchscreen if Analog IO is enabled.

Please contact Thermo Fisher Scientific Sales Service and Customer Support if you need any additional information, see inside cover for contact instructions.

## Testing the Safety Features

The safety features for high temperature protection and low liquid level protection must be checked at regular intervals. We recommend checking at least twice a year or if the bath circulator is moved.

### High temperature protection

Note the pretest position of the arrow on the red dial.

Start the bath and, using a flathead screwdriver, turn the red dial counterclockwise until the bath shuts down and an error message appears.

If the message does not appear or the bath does not shut down, or the message does not clear, have the bath checked by a qualified technician.

Turn the red dial clockwise back to the pretest position.

Allow the bath to cool down before restarting.

To restart the bath press the black reset ring surrounding the red dial. If Auto Restart is enabled the bath restarts, if disabled use the Daily Start Up procedure.

### Low liquid level protection

With the bath on, insert a 1/16" diameter tool, approximately 6" long, into one of the small round openings located at the rear of the reservoir. Slowly push down until the low level error message appears and the bath shuts down. Remove the tool and then clear the error message.

Repeat the procedure using the other round opening.

If the message does not appear or the bath does not shut down, or the message does not clear, have the bath checked by a qualified technician.

### **High liquid level protection**

With the bath on, slowly fill the reservoir above the max level mark until the high level error message appears and the bath shuts down. Drain the reservoir to below the max level mark and then clear the error message.

If the message does not appear or the bath does not shut down, or the message does not clear, have the bath checked by a qualified technician.

## Appendix A AC Serial Communications Protocol

**Note** Currently serial communications is sent using only the connections on the rear of the bath.

**Note** This appendix assumes you have a basic understanding of communications protocols. Information on the NC, STANDARD and NAMUR protocols is available upon request.

**Note** NC protocol is required to use RS485 device addressing.

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format and responses. Bath responses are either the requested data or an error message. The bath response must be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the bath's response. If the command is not understood, the bath responds with an error command. Otherwise, the bath responds with the requested data.

Commands are not case sensitive. Upper or lower case letters may be used. Commands are listed in the Commands Table, error responses are given in the Errors Table, and symbols are shown in the Key Table.

| Key    |   |
|--------|---|
| Symbol | Meaning   |
| [B]    | A binary value 0 or 1 (0 = Off, FALSE or Disable(d); 1 = On, TRUE or Enable(d)).  |
| [CR]   | Carriage return – used as the termination character.                              |
| [U]    | Text representing the units associated with a value.                              |
| [V]    | A value that can be requested in a read command or sent as part of a set command. |

Value: Read commands return analog [V] or bit [B] values or settings, set commands send analog or bit settings. Read commands return values in the same displayed precision. Set command messages missing the space character between the command and the setting are rejected, as the user's intent is unclear.

Units: A read command returning an analog [V] value or setting, will include the units [U] associated with that value or setting. A set command sending an analog value will not include the units. The units returned by the complementary read command are assumed.

Termination character: A carriage return [CR] is used to terminate command and response messages. (Typically the "Enter" key on the keyboard.)

**Note** The inter-character time out (time between transmitted characters) is set to 15 seconds. Exceeding the time out will clear the receiver buffer and require the message to be retransmitted.

**Note** Special characters (backspace, delete, insert, etc.) are not recognized and generate error responses.

| Commands Table  |              |                                   |
|---|--------------|-----------------------------------|
| Commands  |              |                                   |
| All messages from master and slave are terminated with a carriage return [CR] |              |                                   |
| Command Description   | Master Sends | Sample Slave Response (echo off)* |
| Read Control Temperature with RTA   | RT           | [V]C or F                         |
| Read Remote Temperature with RTA  | RT2          | [V]C or F                         |
| Read Internal Temperature with RTA  | RT3          | [V]C or F                         |
| Read Outlet Temperature with RTA  | RT4          | [V]C or F                         |
| Read Set Point  | RS           | [V]C or F                         |
| Read Set Point Preset 1   | RS1          | [V]C or F                         |
| Read Set Point Preset 2   | RS2          | [V]C or F                         |
| Read Set Point Preset 3   | RS3          | [V]C or F                         |
| Read Set Point Preset 4   | RS4          | [V]C or F                         |
| Read Set Point Preset 5   | RS5          | [V]C or F                         |
| Read Internal RTA   | RIRTA        | [V]C or F                         |
| Read Internal RTA Preset 1  | RIRTA1       | [V]C or F                         |
| Read Internal RTA Preset 2  | RIRTA2       | [V]C or F                         |
| Read Internal RTA Preset 3  | RIRTA3       | [V]C or F                         |
| Read Internal RTA Preset 4  | RIRTA4       | [V]C or F                         |
| Read Internal RTA Preset 5  | RIRTA5       | [V]C or F                         |
| Read External RTA   | RERTA        | [V]C or F                         |
| Read External RTA Preset 1  | RERTA1       | [V]C or F                         |
| Read External RTA Preset 2  | RERTA2       | [V]C or F                         |
| Read External RTA Preset 3  | RERTA3       | [V]C or F                         |
| Read External RTA Preset 4  | RERTA4       | [V]C or F                         |
| Read External RTA Preset 5  | RERTA5       | [V]C or F                         |

\* see Possible Response Messages on page A-5.

| Command Description                    | Master Sends | Sample Slave Response (echo off)*   |
|--|--------------|---|
| Read Temperature Units                 | RTU          | [V]C or F   |
| Read Temperature Resolution            | RTR          | [V]   |
| Read High Temperature Limit            | RHT          | [V]C or F   |
| Read High Temperature Fault            | RHTF         | [V]C or F   |
| Read High Temperature Warning          | RHTW         | [V]C or F   |
| Read Low Temperature Warning           | RLTW         | [V]C or F   |
| Read Low Temperature Fault             | RLTF         | [V]C or F   |
| Read Low Temperature Limit             | RLT          | [V]C or F   |
| Read Proportional Heat Band Setting    | RPH          | [V]   |
| Read Integral Heat Band Setting        | RIH          | [V]   |
| Read Derivative Heat Band Setting      | RDH          | [V]   |
| Read Bath On Status                    | RO           | [B]   |
| Read Internal Control on Outlet Sensor | RINTCTL      | [V]<br>V=1 Internal Sensor,<br>V=2 Outlet Sensor  |
| Read External Sensor Enabled           | RE           | [B]   |
| Read Auto Restart Enabled              | RAR          | [B]   |
| Read Pump Speed                        | RPS          | [V] (V is 40 to 100 in %)   |
| Read Fluid Selection                   | RFLUIDTYP    | [V] V is 0=H <sub>2</sub> O, 1=EG/H <sub>2</sub> O,<br>2=PG/H <sub>2</sub> O, 3=SIL 180, 4=SIL200 |
| Read Firmware Version                  | RVER         | [V]   |
| Read Firmware Checksum                 | RSUM         | [V]   |
| Read Unit Fault Status                 | RUFS         | [V1, V2 , V3, V4]<br>See page A-4   |

\* see Possible Response Messages on page A-5.

| Commands Table   |                     |                                  |
|--|---------------------|----------------------------------|
| Commands All messages from master and slave are terminated with a carriage return [CR] |                     |                                  |
| Command Description  | Master Sends        | Sample Slave Response (echo off) |
| Set Set Point  | SS [V]              | OK                               |
| Set Set Point Preset 1   | SS1 [V]             | OK                               |
| Set Set Point Preset 2   | SS2 [V]             | OK                               |
| Set Set Point Preset 3   | SS3 [V]             | OK                               |
| Set Set Point Preset 4   | SS4 [V]             | OK                               |
| Set Set Point Preset 5   | SS5 [V]             | OK                               |
| Set Internal RTA   | SIRTA               | OK                               |
| Set Internal RTA Preset 1  | SIRTA1 [V]          | OK                               |
| Set Internal RTA Preset 2  | SIRTA2 [V]          | OK                               |
| Set Internal RTA Preset 3  | SIRTA3 [V]          | OK                               |
| Set Internal RTA Preset 4  | SIRTA4 [V]          | OK                               |
| Set Internal RTA Preset 5  | SIRTA5 [V]          | OK                               |
| Set External RTA   | SERTA               | OK                               |
| Set External RTA Preset 1  | SERTA1 [V]          | OK                               |
| Set External RTA Preset 2  | SERTA2 [V]          | OK                               |
| Set External RTA Preset 3  | SERTA3 [V]          | OK                               |
| Set External RTA Preset 4  | SERTA4 [V]          | OK                               |
| Set External RTA Preset 5  | SERTA5 [V]          | OK                               |
| Set High Temperature Fault   | SHTF [V]            | OK                               |
| Set High Temperature Warning   | SHTW [V]            | OK                               |
| Set Low Temperature Fault  | SLTF [V]            | OK                               |
| Set Low Temperature Warning  | SLTW [V]            | OK                               |
| Set Proportional Heat Band Setting   | SPH [V]             | OK                               |
| Set Integral Heat Band Setting   | SIH [V]             | OK                               |
| Set Derivative Heat Band Setting   | SDH [V]             | OK                               |
| Set Temperature Resolution   | STR [V]             | OK                               |
| Set Temperature Units  | STU [V] V is C or F | OK                               |
| Set Bath On Status   | SO [B]              | OK                               |
| Set Internal Control on Outlet Sensor<br>V=1 Internal Sensor, V=2 Outlet Sensor        | SINTCTL [V]         | OK                               |
| Set External Sensor Enabled  | SE [B]              | OK                               |
| Set Auto Restart Enabled   | SAR [B]             | OK                               |
| Set Pump Speed   | SPS [V]             | OK                               |
| Set Unit Fault Status (clear errors)   | SUFS                | OK                               |

**RUFS (Read Unit Fault Status)**

This command returns 4 hex values (e.g., 00 00 00 B8) – bath on with pump, compressor and heater running.

| hex | B7 | B6 | B5 | B4 | B3 | B2 | B1 | B0 |
|-----|----|----|----|----|----|----|----|----|
| 1   | 0  | 0  | 0  | 0  | 0  | 0  | 0  | 1  |
| 2   | 0  | 0  | 0  | 0  | 0  | 0  | 1  | 0  |
| 4   | 0  | 0  | 0  | 0  | 0  | 1  | 0  | 0  |
| 8   | 0  | 0  | 0  | 0  | 1  | 0  | 0  | 0  |
| 10  | 0  | 0  | 0  | 1  | 0  | 0  | 0  | 0  |
| 20  | 0  | 0  | 1  | 0  | 0  | 0  | 0  | 0  |
| 40  | 0  | 1  | 0  | 0  | 0  | 0  | 0  | 0  |
| 80  | 1  | 0  | 0  | 0  | 0  | 0  | 0  | 0  |

| Value | Description of bits   |
|-------|---|
| V1    | B7 rtd1 open<br>B6 rtd1 shorted<br>B5 rtd2 open<br>B4 rtd2 shorted<br>B3 rtd3 open<br>B2 rtd3 shorted<br>B1 remote rtd open<br>B0 remote rtd shorted                          |
| V2    | B7 high temp fixed fault<br>B6 low temp fixed fault<br>B5 high temp fault<br>B4 low temp fault<br>B3 high ra temp fixed fault<br>B2 HTC fault<br>B1 HPC fault<br>B0 LLC fault |
| V3    | B7 low level warning<br>B6 high temp warning<br>B5 low temp warning<br>B4 Auto Refill Timeout warning<br>B3 min step out of range warning<br>B2 - B0 0 (reserved)             |
| V4    | B7 Bath on<br>B6 Bath Faulted<br>B5 Pump Relay on<br>B4 Compressor Relay on<br>B3 Heating<br>B2 - B0 0 (reserved)   |

Refer to Key table on page 1 for explanation of symbols and their meanings.

Examples:

Read Temperature:

Host

|         |   |           |
|---------|---|-----------|
| R       | T | <b>CR</b> |
| Command |   | [CR]      |

Bath:

|     |   |   |   |     |           |
|-----|---|---|---|-----|-----------|
| 2   | 0 | . | 0 | C   | <b>CR</b> |
| [V] |   |   |   | [U] | [CR]      |

Set Set point:

Host

|         |   |  |     |      |           |
|---------|---|--|-----|------|-----------|
| S       | S |  | 2   | 0    | <b>CR</b> |
| Command |   |  | [V] | [CR] |           |

Bath:

|                  |   |           |
|------------------|---|-----------|
| O                | K | <b>CR</b> |
| Command Accepted |   | [CR]      |

Read Temperature 2:

Host:

|       |   |   |           |   |   |      |
|-------|---|---|-----------|---|---|------|
| R     | T | 2 | <b>CR</b> |   |   |      |
| Bath: | 2 | 0 | .         | 0 | C | [CR] |

Set Set point to -22°C when minimum allowed is -20°C: Minimum allowed is [VMIN]


Host:

|       |   |   |   |   |   |           |  |   |   |   |  |   |   |  |   |   |   |   |   |           |
|-------|---|---|---|---|---|-----------|--|---|---|---|--|---|---|--|---|---|---|---|---|-----------|
| S     | S |   | - | 2 | 2 | <b>CR</b> |  |   |   |   |  |   |   |  |   |   |   |   |   |           |
| Bath: | F | A | I | L |   | -         |  | O | u | t |  | o | f |  | R | a | n | g | e | <b>CR</b> |

| Possible Response Messages |   |
|----------------------------|---|
| "OK"                       | Command accepted OK   |
| "FAIL"                     | Service: Calibration failed                                 |
| "FAIL - Invalid Setting"   | For a small set of discrete values like 1/0 (on/off); C/F   |
| "FAIL - Out of Range"      | For continuous settings such as a set point                 |
| "FAIL - Password Required" | Service: Must enter a password to use (e.g., set fan speed) |
| "?"                        | Unknown command or not implemented                          |


# B

## Appendix B Analog I/O


Install your analog input/output device to the 15-pin female connector on the rear of the bath. It is activated using the touchscreen.

**Note** You cannot save changes using the touchscreen when Analog I/O is enabled.

Use the **SetPoint IN** and **Temp OUT** screens to select the analog input and DAC output scaling to meet the application's needs.

The VersaCool supports three standard analog interface types and defaults using the following scaling:  
**mV** 10mv/°C where 0v = 0°C, e.g., 50mV = +5°C, -100mV = -10°C  
**Volt** where 0v – 10v = the bath's operating range (-20°C to +150°C)  
**mA** where 4ma – 20ma = the bath's operating range (-20°C to +150°C)


Use **Custom** to change the factory defaults and enter your own ranges for any of the three options. The factory defaults are:

**mA** where 4ma – 20ma = the bath's operating range (-20°C to +150°C)

**Volt** where 0v – 10v = the bath's operating range (-20°C to +150°C)

**mV** 10mv/°C where 0V = 0°C, e.g., 50mV = +5°C, -100mV = -10°C

Use the **Relay Status** screen to select the desired error(s).

**Relay 1**


- Any Fault
- Low Level Fault
- Low Temperature
- High Temperature
- Low Level Warning
- Pump/Unit Shutdown
- RA Shutdown
- Limit Fault
- Sensor Fault

**Relay 2**

- Any Fault
- Low Level Fault
- Low Temperature
- High Temperature
- Low Level Warning
- Pump/Unit Shutdown
- RA Shutdown
- Limit Fault
- Sensor Fault

**Never apply line voltage to any of the connections.**

| Analog IO Connector Pinout |   |       |   |
|----------------------------|---|-------|---|
| Pin                        | Name  | Notes | Definition  |
| 1                          | Digital ground  |       | Common round connection for pins 12, 13 and 14  |
| 2                          | Not used  |       |   |
| 3                          | Low level   | 1     | Dry Relay Contact: Reference to pin 11.<br>Closes if either level switch is in the "low" position for more than 1 second.                 |
| 4                          | Configurable relay 2  | 1     | <u>Dry Relay Contact</u> : Reference to pin 11.<br>Closes when any configured fault or warning occurs.                                    |
| 5                          | Pump on   | 1     | <u>Dry Relay Contact</u> : Reference to pin 11.<br>Closes when pump is turned on. opens when pump is off                                  |
| 6                          | Analog ground   |       | Common for analog signals (pins 2, 7 and 15)  |
| 7                          | Reservoir temp out <b>or</b> remote sensor temp if remote temp sensor enabled | 2     | Analog Voltage Output: Reference to pin 6.  |
| 8                          | High level  | 1     | <u>Dry Relay Contact</u> : Reference to pin 11<br>Closes if level switch is in the high position for approximately 1 second or more.      |
| 9                          | Configurable relay 1 (normally open)  | 1     | <u>Dry Relay Contact</u> : Reference to pin 11<br>Closes when any of the configured faults occur.   |
| 10                         | Configurable relay 1 (normally closed)  | 1     | <u>Dry Relay Contact</u> : Reference to pin 11<br>Complement of pin 9 (open when pin 9 is closed).  |
| 11                         | Relay common  |       | Common for all relay contacts (pins 3, 4, 5, 8, 9, 10).   |
| 12                         | Not used  |       |   |
| 13                         | Remote set point enabled  | 3     | Connect to pin 1 to allow the set point to be changed remotely through pin 15 Remote set point.   |
| 14                         | Remote start  | 3     | Connect to pin 1 to turn bath on. Disconnect to turn bath off.<br>Note: Pins 1 and 12 must be connected to allow operation from this pin. |
| 15                         | Remote set point  | 2, 4  | Analog Voltage Input: Reference to pin 6.   |

Note 1: All relay contacts (except for Pin 10) are normally OPEN when power is off. Pin 10 contacts are normally CLOSED when power is off. Relay contacts are rated: 24V AC/DC, 2A, <= 0.08 Ohm maximum each or 5A total for all relays combined, 1mA minimum, switching capacity: 48VA/48W (resistive load only).


Note 2: Default = 0-10VDC. These ranges are set by the user.

Note 3: Connect to digital ground (pin 1) using a low resistance connection (gold contact relay).

Note 4: Remote set point must be enabled, pin 13

When making your connection to the VersaCool Analog I/O connector, in order to comply with the EMC directive:

- Use a shielded I/O cable
- Connect the remote end of the cable shield to earth ground.
- Connect cable shield to VersaCool end connector.

# C


## Appendix C Mobile Communications

The VersaCool Mobile App can use any iOS™ or Android™ smart phone or tablet to control and monitor up to seven VersaCool baths. Mobile communications is enabled when the circuit protector on the rear of the bath is in the **I** position.

The Mobile App allows you to turn multiple VersaCool's on/off, view their fluid temperature at the bath's selected probe, enter custom set points, select preset set points, view alarm statuses, turn alarm audio on/off, and monitor a chart for up to seven baths.


After start up touch the red dot and the mobile app will locate any powered VersaCool's within range. A spinning circle appears when the search is underway.


Once activated the mobile app will search for and then list all found Bluetooth Devices. The list includes the VersaCool's Bluetooth address.

**Note** A VersaCool cannot be selected by more than one mobile device at a time.

A message appears if no bath is found.

If you are not sure which baths are on the list refer to the VersaCool's **Configuration** display to identify its **Bluetooth MAC** address.


Touch the desired bath's name to display a drop down showing the bath's address, pump status, last recorded reservoir fluid temperature, current set point and current time and date, bath reservoir fluid temperature, set point and alarm history.

This display will also indicate any alarms. Touch the alarm to remove it from the display.

**Note** You cannot clear the alarm using the device. For safety, the alarm can only be cleared at the bath.

Touch  to start the bath.

The red bar will turn green and now display , see below .


Touch **Enter custom SP** to bring up the keypad used to change the set point. Touch **Choose SP1 - SP5** to select a preset set point. SP 1– 5 are the same set points stored on the displayed VersaCool.

**Note** These preset set points can only be selected, not changed, using the mobile app.


Touching  stops the bath.

Touch **< VersaCool Mobile** at the top-left of the display to return to the list of available baths. Use this list to activate/deactivate any of the displayed baths. See next page.


Touch  to bring up a list of all the baths currently connected to your mobile device, their reservoir fluid temperatures and their statuses.

Touch  to bring up a list of Icon Descriptions.

Touch  to bring up a graphical display of the bath's performance.

The chart feature displays up to 30 minutes of past performance. 0 minutes is on the right, 30 minutes is on the left.

Each bath is identified with a different color.

Touch  to toggle between displaying that bath or not.

**Note** Any break in the chart indicates where connectivity was lost.


Touch  to pause charting.


Touch  to send the displayed image as a screen shot over e-mail or as a text message (as long as the device is connected to your network).

Rotate your mobile device to bring up a full screen display, see below.

Pinch and drag any portion of the chart to zoom in.

**Note** The chart data is for display purposes only, it is not permanently stored on the mobile device. If you wish to access historical data you must download the data directly from the VersaCool (e.g. USB stick, communication port).


Touch  to display the chart's temperature settings.


You can select the desired temperature scale, Fahrenheit or Celsius.

Touch the displayed min or max temp value to bring up the keypad you use to make changes.


If desired, touch and hold the bath's name to bring up the swatch used to change its color. Touch the desired color.

Cancel or Ok your selection.

Touch **⋮** on the right of the desired bath to display its Instrument Settings.

Use this display to rename the bath, enable/disable its audible alarms, connect/disconnect the bath and/or remove it from the active list.

**Note** Changing the name only affects how it appears on your mobile device. It does not change the name at the bath itself.

## **GNU LESSER GENERAL PUBLIC LICENSE**

Version 2.1, February 1999

Copyright (C) 1991, 1999 Free Software Foundation, Inc.

51 Franklin Street, Fifth Floor

Boston, MA 02110-1301 USA

Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

BECAUSE THE LIBRARY IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE LIBRARY, TO THE EXTENT PERMITTED BY APPLICABLE LAW.

The Open Source code used in the development of this product will be provided for a period of 3-years from date of purchase to the original purchaser provided they supply Thermo Fisher Scientific with the applicable serial number for the product the Open Source code is being requested for.

[This is the first released version of the Lesser GPL. It also counts as the successor of the GNU Library Public License, version 2, hence the version number 2.1.]

### Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public Licenses are intended to guarantee your freedom to share and change free software--to make sure the software is free for all its users.

This license, the Lesser General Public License, applies to some specially designated software packages--typically libraries--of the Free Software Foundation and other authors who decide to use it. You can use it too, but we suggest you first think carefully about whether this license or the ordinary General Public License is the better strategy to use in any particular case, based on the explanations below.

When we speak of free software, we are referring to freedom of use, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish); that you receive source code or can get it if you want it; that you can change the software and use pieces of it in new free programs; and that you are informed that you can do these things.

To protect your rights, we need to make restrictions that forbid distributors to deny you these rights or to ask you to surrender these rights. These restrictions translate to certain responsibilities for you if you distribute copies of the library or if you modify it.

For example, if you distribute copies of the library, whether gratis or for a fee, you must give the recipients all the rights that we gave you. You must make sure that they, too, receive or can get the source code. If you link other code with the library, you must provide complete object files to the recipients, so that they can relink them with the library after making changes to the library and recompiling it. And you must show them these terms so they know their rights.

We protect your rights with a two-step method: (1) we copyright the library, and (2) we offer you this license, which gives you legal permission to copy, distribute and/or modify the library.

To protect each distributor, we want to make it very clear that there is no warranty for the free library. Also, if the library is modified by someone else and passed on, the recipients should know that what they have is not the original version, so that the original author's reputation will not be affected by problems that might be introduced by others.

Finally, software patents pose a constant threat to the existence of any free program. We wish to make sure that a company cannot effectively restrict the users of a free program by obtaining a restrictive license from a patent holder. Therefore, we insist that any patent license obtained for a version of the library must be consistent with the full freedom of use specified in this license.

Most GNU software, including some libraries, is covered by the ordinary GNU General Public License. This license, the GNU Lesser General Public License, applies to certain designated libraries, and is quite different from the ordinary General Public License. We use this license for certain libraries in order to permit linking those libraries into non-free programs.

When a program is linked with a library, whether statically or using a shared library, the combination of the two is legally speaking a combined work, a derivative of the original library. The ordinary General Public License therefore permits such linking only if the entire combination fits its criteria of freedom. The Lesser General Public License permits more lax criteria for linking other code with the library.

We call this license the "Lesser" General Public License because it does Less to protect the user's freedom than the ordinary General Public License. It also provides other free software developers Less of an advantage over competing non-free programs. These disadvantages are the reason we use the ordinary General Public License for many libraries. However, the Lesser license provides advantages in certain special circumstances.

For example, on rare occasions, there may be a special need to encourage the widest possible use of a certain library, so that it becomes a de-facto standard. To achieve this, non-free programs must be allowed to use the library. A more frequent case is that a free library does the same job as widely used non-free libraries. In this case, there is little to gain by limiting the free library to free software only, so we use the Lesser General Public License.

In other cases, permission to use a particular library in non-free programs enables a greater number of people to use a large body of free software. For example, permission to use the GNU C Library in non-free programs enables many more people to use the whole GNU operating system, as well as its variant, the GNU/Linux operating system.

Although the Lesser General Public License is Less protective of the users' freedom, it does ensure that the user of a program that is linked with the Library has the freedom and the wherewithal to run that program using a modified version of the Library.

The precise terms and conditions for copying, distribution and modification follow. Pay close attention to the difference between a "work based on the library" and a "work that uses the library". The former contains code derived from the library, whereas the latter must be combined with the library in order to run.

## GNU LESSER GENERAL PUBLIC LICENSE TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

0. This License Agreement applies to any software library or other program which contains a notice placed by the copyright holder or other authorized party saying it may be distributed under the terms of this Lesser General Public License (also called "this License"). Each licensee is addressed as "you".

A "library" means a collection of software functions and/or data prepared so as to be conveniently linked with application programs (which use some of those functions and data) to form executables.

The "Library", below, refers to any such software library or work which has been distributed under these terms. A "work based on the Library" means either the Library or any derivative work under copyright law: that is to say, a work containing the Library or a portion of it, either verbatim or with modifications and/or translated straightforwardly into another language. (Hereinafter, translation is included without limitation in the term "modification".)

"Source code" for a work means the preferred form of the work for making modifications to it. For a library, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the library.

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running a program using the Library is not restricted, and output from such a program is covered only if its contents constitute a work based on the Library (independent of the use of the Library in a tool for writing it). Whether that is true depends on what the Library does and what the program that uses the Library does.

1. You may copy and distribute verbatim copies of the Library's complete source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and distribute a copy of this License along with the Library.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

2. You may modify your copy or copies of the Library or any portion of it, thus forming a work based on the Library, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:

- a) The modified work must itself be a software library.
- b) You must cause the files modified to carry prominent notices stating that you changed the files and the date of any change.
- c) You must cause the whole of the work to be licensed at no charge to all third parties under the terms of this License.
- d) If a facility in the modified Library refers to a function or a table of data to be supplied by an application program that uses the facility, other than as an argument passed when the facility is invoked, then you must make a good faith effort to ensure that, in the event an application does not supply such function or table, the facility still operates, and performs whatever part of its purpose remains meaningful.

(For example, a function in a library to compute square roots has a purpose that is entirely well-defined independent of the application. Therefore, Subsection 2d requires that any application-supplied function or table used by this function must be optional: if the application does not supply it, the square root function must still compute square roots.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Library, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Library, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Library.

In addition, mere aggregation of another work not based on the Library with the Library (or with a work based on the Library) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

3. You may opt to apply the terms of the ordinary GNU General Public License instead of this License to a given copy of the Library. To do this, you must alter all the notices that refer to this License, so that they refer to the ordinary GNU General Public License, version 2, instead of to this License. (If a newer version than version 2 of the ordinary GNU General Public License has appeared, then you can specify that version instead if you wish.) Do not make any other change in these notices.

Once this change is made in a given copy, it is irreversible for that copy, so the ordinary GNU General Public License applies to all subsequent copies and derivative works made from that copy. This option is useful when you wish to copy part of the code of the Library into a program that is not a library.

4. You may copy and distribute the Library (or a portion or derivative of it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange.

If distribution of object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place satisfies the requirement to distribute the source code, even though third parties are not compelled to copy the source along with the object code.

5. A program that contains no derivative of any portion of the Library, but is designed to work with the Library by being compiled or linked with it, is called a "work that uses the Library". Such a work, in isolation, is not a derivative work of the Library, and therefore falls outside the scope of this License.

However, linking a "work that uses the Library" with the Library creates an executable that is a derivative of the Library (because it contains portions of the Library), rather than a "work that uses the library". The executable is therefore covered by this License. Section 6 states terms for distribution of such executables.

When a "work that uses the Library" uses material from a header file that is part of the Library, the object code for the work may be a derivative work of the Library even though the source code is not. Whether this is true is especially significant if the work can be linked without the Library, or if the work is itself a library. The threshold for this to be true is not precisely defined by law.

If such an object file uses only numerical parameters, data structure layouts and accessors, and small macros and small inline functions (ten lines or less in length), then the use of the object file is unrestricted, regardless of whether it is legally a derivative work. (Executables containing this object code plus portions of the Library will still fall under Section 6.)

Otherwise, if the work is a derivative of the Library, you may distribute the object code for the work under the terms of Section 6. Any executables containing that work also fall under Section 6, whether or not they are linked directly with the Library itself.

6. As an exception to the Sections above, you may also combine or link a "work that uses the Library" with the Library to produce a work containing portions of the Library, and distribute that work under terms of your choice, provided that the terms permit modification of the work for the customer's own use and reverse engineering for debugging such modifications.

You must give prominent notice with each copy of the work that the Library is used in it and that the Library and its use are covered by this License. You must supply a copy of this License. If the work during execution displays copyright notices, you must include the copyright notice for the Library among them, as well as a reference directing the user to the copy of this License. Also, you must do one of these things:

a) Accompany the work with the complete corresponding machine-readable source code for the Library including whatever changes were used in the work (which must be distributed under Sections 1 and 2 above); and, if the work is an executable linked with the Library, with the complete machine-readable "work that uses the Library", as object code and/or source code, so that the user can modify the Library and then relink to produce a modified executable containing the modified Library. (It is understood that the user who changes the contents of definitions files in the Library will not necessarily be able to recompile the application to use the modified definitions.)

b) Use a suitable shared library mechanism for linking with the Library. A suitable mechanism is one that (1) uses at run time a copy of the library already present on the user's computer system, rather than copying library functions into the executable, and (2) will operate properly with a modified version of the library, if the user installs one, as long as the modified version is interface-compatible with the version that the work was made with.

c) Accompany the work with a written offer, valid for at least three years, to give the same user the materials specified in Subsection 6a, above, for a charge no more than the cost of performing this distribution.

d) If distribution of the work is made by offering access to copy from a designated place, offer equivalent access to copy the above specified materials from the same place.

e) Verify that the user has already received a copy of these materials or that you have already sent this user a copy.

For an executable, the required form of the "work that uses the Library" must include any data and utility programs needed for reproducing the executable from it. However, as a special exception, the materials to be distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable.

It may happen that this requirement contradicts the license restrictions of other proprietary libraries that do not normally accompany the operating system. Such a contradiction means you cannot use both them and the Library together in an executable that you distribute.

7. You may place library facilities that are a work based on the Library side-by-side in a single library together with other library facilities not covered by this License, and distribute such a combined library, provided that the separate distribution of the work based on the Library and of the other library facilities is otherwise permitted, and provided that you do these two things:

a) Accompany the combined library with a copy of the same work based on the Library, uncombined with any other library facilities. This must be distributed under the terms of the Sections above.

b) Give prominent notice with the combined library of the fact that part of it is a work based on the Library, and explaining where to find the accompanying uncombined form of the same work.

8. You may not copy, modify, sublicense, link with, or distribute the Library except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense, link with, or distribute the Library is void, and will automatically terminate your rights under this License. However, parties who have received copies, or rights, from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

9. You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Library or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Library (or any work based on the Library), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Library or works based on it.

10. Each time you redistribute the Library (or any work based on the Library), the recipient automatically receives a license from the original licensor to copy, distribute, link with or modify the Library subject to these terms and conditions. You may not impose any further restrictions on the recipients' exercise of the rights granted herein. You are not responsible for enforcing compliance by third parties with this License.

11. If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Library at all. For example, if a patent license would not permit royalty-free redistribution of the Library by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Library.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply, and the section as a whole is intended to apply in other circumstances.

It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

12. If the distribution and/or use of the Library is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Library under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.

13. The Free Software Foundation may publish revised and/or new versions of the Lesser General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns.

Each version is given a distinguishing version number. If the Library specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by

the Free Software Foundation. If the Library does not specify a license version number, you may choose any version ever published by the Free Software Foundation.

14. If you wish to incorporate parts of the Library into other free programs whose distribution conditions are incompatible with these, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

#### NO WARRANTY

15. BECAUSE THE LIBRARY IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE LIBRARY, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE LIBRARY "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE LIBRARY IS WITH YOU. SHOULD THE LIBRARY PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

16. IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE LIBRARY AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE LIBRARY (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE LIBRARY TO OPERATE WITH ANY OTHER SOFTWARE), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

#### END OF TERMS AND CONDITIONS

##### How to Apply These Terms to Your New Libraries

If you develop a new library, and you want it to be of the greatest possible use to the public, we recommend making it free software that everyone can redistribute and change. You can do so by permitting redistribution under these terms (or, alternatively, under the terms of the ordinary General Public License).

To apply these terms, attach the following notices to the library. It is safest to attach them to the start of each source file to most effectively convey the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

<one line to give the library's name and a brief idea of what it does.> Copyright (C) <year> <name of author>

This library is free software; you can redistribute it and/or modify it under the terms of the GNU Lesser General Public License as published by the Free Software Foundation; either version 2.1 of the License, or (at your option) any later version.

This library is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU Lesser General Public License for more details.

You should have received a copy of the GNU Lesser General Public License along with this library; if not, write to the Free Software Foundation, Inc., 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301 USA

Also add information on how to contact you by electronic and paper mail.

You should also get your employer (if you work as a programmer) or your school, if any, to sign a "copyright disclaimer" for the library, if necessary. Here is a sample; alter the names:

Yoyodyne, Inc., hereby disclaims all copyright interest in the library `Frob' (a library for tweaking knobs) written by James Random Hacker.

<signature of Ty Coon>, 1 April 1990

Ty Coon, President of Vice

That's all there is to it!

Nokia Qt LGPL Exception version 1.1

As an additional permission to the GNU Lesser General Public License version 2.1, the object code form of a "work that uses the Library" may incorporate material from a header file that is part of the Library. You may distribute such object code under terms of your choice, provided that:

- (i) the header files of the Library have not been modified; and
- (ii) the incorporated material is limited to numerical parameters, data structure layouts, accessors, macros, inline functions and templates; and
- (iii) you comply with the terms of Section 6 of the GNU Lesser General Public License version 2.1.

Moreover, you may apply this exception to a modified version of the Library, provided that such modification does not involve copying material from the Library into the modified Library's header files unless such material is limited to (i) numerical parameters; (ii) data structure layouts; (iii) accessors; and (iv) small macros, templates and inline functions of five lines or less in length.

Furthermore, you are not required to apply this additional permission to a modified version of the Library.

# Warranty

Thermo Fisher Scientific warrants electronic components for 36 months, 24 months on all other components, from date of shipment the Thermo Scientific VersaCool Refrigerated Bath according to the following terms.

Any part of the bath manufactured or supplied by Thermo Fisher Scientific and found in the reasonable judgment of Thermo Fisher to be defective in material or workmanship will be repaired at an authorized Thermo Fisher Repair Depot without charge for parts or labor. The bath, including any defective part must be returned to an authorized Thermo Fisher Repair Depot within the warranty period. The expense of returning the bath to the authorized Thermo Fisher Repair Depot for warranty service will be paid for by the buyer. Our responsibility in respect to warranty claims is limited to performing the required repairs or replacements, and no claim of breach of warranty shall be cause for cancellation or rescission of the contract of sales of any bath. With respect to baths that qualify for field service repairs, Thermo Fisher Scientific's responsibility is limited to the component parts necessary for the repair and the labor that is required on site to perform the repair. Any travel labor or mileage charges are the financial responsibility of the buyer.

The buyer shall be responsible for any evaluation or warranty service call (including labor charges) if no defects are found with the Thermo Scientific product.

This warranty does not cover any bath that has been subject to misuse, neglect, or accident. This warranty does not apply to any damage to the bath that is the result of improper installation or maintenance, or to any bath that has been operated or maintained in any way contrary to the operating or maintenance instructions specified in this Instruction and Operation Manual. This warranty does not cover any bath that has been altered or modified so as to change its intended use.

In addition, this warranty does not extend to repairs made by the use of parts, accessories, or fluids which are either incompatible with the bath or adversely affect its operation, performance, or durability.

Thermo Fisher Scientific reserves the right to change or improve the design of any bath without assuming any obligation to modify any bath previously manufactured.

THE FOREGOING EXPRESS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

OUR OBLIGATION UNDER THIS WARRANTY IS STRICTLY AND EXCLUSIVELY LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE COMPONENT PARTS AND Thermo Fisher Scientific DOES NOT ASSUME OR AUTHORIZE ANYONE TO ASSUME FOR IT ANY OTHER OBLIGATION.

Thermo Fisher Scientific ASSUMES NO RESPONSIBILITY FOR INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES INCLUDING, BUT NOT LIMITED TO LOSS OR DAMAGE TO PROPERTY, LOSS OF PROFITS OR REVENUE, LOSS OF THE BATH, LOSS OF TIME, OR INCONVENIENCE.

This warranty applies to baths sold by Thermo Fisher Scientific. (Refer to the warranty for baths sold by the affiliated marketing company of Thermo Fisher Scientific for any additional terms.) This warranty and all matters arising pursuant to it shall be governed by the law of the State of New Hampshire, United States. All legal actions brought in relation hereto shall be filed in the appropriate state or federal courts in New Hampshire, unless waived by Thermo Fisher Scientific.

[thermoscientific.com](http://thermoscientific.com)

© 2016 Thermo Fisher Scientific Inc. All rights reserved. iOS is a registered trademark of Cisco. Android is a registered trademark of Google Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Thermo Fisher Scientific**  
81 Wyman Street  
Waltham, MA 02451

**thermo**  
scientific

A Thermo Fisher Scientific Brand