

# Powerful at Every Turn


Thermo Scientific™ General Purpose Pro Centrifuges



## New Touchscreen Technology

**Easy menu access:** no more scrolling

**Pre-temp:** centrifuge can be “at temp” before run



## Home Screen

**Centrifuge Health:** preventative maintenance alerts and cues for less down time

**Help menu:** virtual manual provides quick solutions

**Pulse:** customizable one-touch quick spins



## New Auto-Lock™ II Rotor Exchange



Quick switch between applications



Seconds to safely switch rotors



New soft unlock feature



No tools required

## FiberLite™ Carbon Fiber Rotors

- Lightweight carbon fiber
- Outstanding durability
- Corrosion and fatigue resistant
- Backed by a 15-year limited warranty



## ClickSeal™ Biocontainment Lids



Glove friendly



One-handed closure



Biocontainment certifications



Designed for right- or left-hand operation



## Capacity and Performance

- Spin up to 40 x 50mL conical tubes or 196 blood sample tubes
- Compact footprint
- Improved sample safety
- Greater sample reproducibility

## Improved Ergonomics



Loading height decreased by 1.2cm (0.5 in.)



Less effort required to close lid



Lid grabbing height decreased by 5.1cm (2 in.)