

FAT Binder Matrix

Summary of factory acceptance test documentation for cell therapy applications

Samenvatting en overzicht van 'Factory acceptance test' documentatie voor celtherapie applicaties

Binder Element	CO ₂	BSC	Cent	CRF	Description	Value to GMP customer
Product support						
Design qualification DQ (Technical Data Sheet, drawings)*					Basic spec descriptions of the product an may include some summary drawings	Helps customers choose product and also define installation requirements
Recommended customer maintenance schedule					Vendor-supplied customer-based maintenance & schedule	Helpful to create time based plan for maintenance and audits
Consumables and customer replaceable parts					Replaceable parts based upon usage or wear	Keeps system up and running with proper attention and scheduling
Recommended cleaning procedures					Effective and typical procedures to keep product clean	Avoids any possibility of contamination
Regulatory and Compliance						
Declaration of conformity with EU laws (CE)					CE certificate proof	Attest to the quality level and specification achievements
CSA Certificate					CSA certificate proof	Attest to the quality level and specification achievements
UL Certificate					UL certificate proof	Attest to the quality level and specification achievements
Certificate of Environmental Status					Regulatory certificates for various sustainability items	Shows compliance to the very dynamic environmental regulations
Quality Management System [QMS] Certificate					Factory/site based Quality system description	Provides evidence of adequate quality procedures and compliance
REACH Statement					Regulatory certificates for various sustainability items	Shows compliance to the very dynamic environmental regulations
DIN EN12469 TUV Certification					Regulatory compliance certificate	Well established clean room and airflow standard
DIN ISO 14644-1 Cleanroom compatibility certificate					Regulatory compliance certificate	Well established clean room and airflow standard

* These documents are available from the factory prior to purchase upon request (in addition to CAD format drawings)

Uw leverancier in NL:



Poly Temp Scientific BV
De Marne 211 Bolsward
info@polytemp.nl
0515-575105



Binder Element	CO ₂	BSC	Cent	CRF	Description	Value to GMP customer
Traceability/Validation						
Certificate of Conformance					Production certificate	Proves product made according to specifications
Factory Acceptance Test Report [FAT]					End of Line test summary with checklist of performance	Important to verify product has passed test at Factory level before shipment
Key Material of Construction Specifications					Technical specifications	Proof of specs such materials, surface finish , etc
Steris VHP compatibility certificate					Test certificate for Vaporized Hydrogen Peroxide cleaning	Defines ability to undergo VHP cleaning
Particle Emission Certificate					Test certificate defining particle emissions of product	Enables customer to comply with their clean room particle levels
Porton Down Certificate					Test certificate for centrifuge rotor biocontainment	Proves bio-hazard level containment
HEPA Filter Certificate					Certificate of filter performance	Required to show proper clean room filtration
Temperature Sensor Calibration document					Certificate from supplier on calibration	Proof of calibration
RH Sensor Calibration Document					Certificate from supplier on calibration	Proof of calibration
CO ₂ Sensor Calibration document					Certificate from supplier on calibration	Proof of calibration
CO ₂ Valve certificate of compliance					Certificate from supplier on calibration	Proof of calibration
Certificate of electro-polishing					Material Certificate of surface finish	Proof of cleanability levels and porosity

Color Denotes Serialized Documents

CTS Series

CO₂ = Thermo Scientific™ Forma™ Steri-cult™ CO₂ Incubator

BSC = Thermo Scientific™ Herasafe™ 2030i Biological Safety Cabinet

CENT= Thermo Scientific™ General Purpose Pro Centrifuge

CRF = Thermo Scientific™ CryoMed™ Controller Rate Freezer

Find out more at thermofisher.com/ctslabequipment

Uw leverancier in NL:



Poly Temp Scientific BV
De Marne 211 Bolsward
info@polytemp.nl
0515-575105

ThermoFisher
SCIENTIFIC

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer application. © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. BRFAT 0820